





FUSIONPlex™ for Ion Torrent—Workflow Overview

Follow the instructions and incubations for each step. All mixing steps should be performed on ice. Pipette up and down 8 times or vortex to mix after re-suspending each lysosphere and spin down prior to incubations and transfers. For incubations, use a lid heated $\geq 100^{\circ}\text{C}$, except where specified otherwise.

Random Priming	Add 20–250 ng RNA to Random Priming reagent for a total volume of 20 μL . Incubate as indicated.	Step	Incubation temperature	Incubation time
		1	65°C	5 min
		2	4°C	Hold
First Strand cDNA Synthesis	Transfer 20 μL Random Priming mixture to First Strand cDNA Synthesis reagent. Incubate as indicated. Transfer 1 μL First Strand cDNA Synthesis mixture to 9 μL water for PreSeq RNA QC Assay.	Step	Incubation temperature	Incubation time
		1	25°C	10 min
		2	42°C	30 min
		3	80°C	20 min
Second Strand cDNA Synthesis	Add 21 μL ultrapure water to 19 μL First Strand cDNA Synthesis mixture. Add 21 μL ultrapure water to 19 μL of First Strand cDNA Synthesis mixture and then transfer 40 μL to Second Strand cDNA Synthesis reagent. Incubate as indicated.  Optional stopping point before purification. Store at -10°C to -30°C .	Step	Incubation temperature	Incubation time
		1	16°C	60 min
		2	75°C	20 min
		3	4°C	Hold
PreSeq RNA QC Assay	Build qPCR reaction and incubate as per your master mix-specific instructions. iTag SYBR Green Supermix (Bio-Rad) – 5 μL 10X VCP Primer Mix – 1 μL Diluted cDNA sample or NTC – 4 μL Incubate as indicated.	Incubation temperature	Incubation time	# of cycles
		95°C	20 [20*] sec	1
		95°C	3 [15*] sec	35
		60°C	30 [60*] sec	
		60–95°C	0.5°C/sec increment	1
* Times in [] are for standard cycling				
End Repair	Transfer 40 μL Second Strand cDNA Synthesis mixture to End Repair reagent. Incubate as indicated with <u>unheated lid</u>	Step	Incubation temperature	Incubation time
		1	25°C	30 min
		2	4°C	Hold
AMPure™ XP clean-up (100 μL). Elute in 20 μL .				
Ligation Step 1	Transfer 20 μL End Repair mixture to Ligation Step 1 reagent. Incubate as indicated.	Step	Incubation temperature	Incubation time
		1	37°C	15 min
		2	4°C	Hold
AMPure™ XP clean-up (50 μL). Elute in 42 μL .				
MBC Adapters	Transfer 40 μL Ligation Step 1 mixture to the MBC adapters.			
Ligation Step 2	Transfer the entire volume of the MBC adapters to Ligation Step 2 reagent. Incubate as indicated with <u>unheated lid</u>.  Optional stopping point before purification. Store at -10°C to -30°C .	Step	Incubation temperature	Incubation time
		1	22°C	5 min
		2	4°C	Hold
		Ligation clean-up beads (50 μL). Elute in 18 μL 5mM NaOH, 75°C 10 min.		

First PCR	Add 2 μ L GSP1 primers to the First PCR reagent. Then add 18 μ L of purified Adapter Ligation mixture and mix. Incubate as indicated.  Optional stopping point after this step. Store purified product (eluted from beads) at -10°C to -30°C .	Incubation temperature	Incubation time	# of cycles
		95 $^{\circ}\text{C}$	3 min	1
		95 $^{\circ}\text{C}$	30 sec	Varies*
		Varies*	10 sec	
		Varies*	Varies*	
		72 $^{\circ}\text{C}$	3 min	1
		4 $^{\circ}\text{C}$	Hold	1
		AMPure™ XP clean-up (24 μ L). Elute in 20 μ L.		

Second PCR	Add 2 μ L GSP2 primers to Second PCR reagent. Then add 18 μ L of purified First PCR mixture and mix. Incubate as indicated.  Optional stopping point after this step. Store purified product (eluted from beads) at -10°C to -30°C .	Incubation temperature	Incubation time	# of cycles
		95 $^{\circ}\text{C}$	3 min	1
		95 $^{\circ}\text{C}$	30 sec	Varies*
		Varies*	10 sec	
		Varies*	Varies*	
		72 $^{\circ}\text{C}$	3 min	1
		4 $^{\circ}\text{C}$	Hold	1
		AMPure™ XP clean-up (24 μ L). Elute in 20 μ L.		

* Refer to Product Insert for panel-specific parameters

Proceed with Ion Torrent™ Library Quantification Assay and Sequence

Appendix: AMPure® XP bead purification†

1. Add AMPure XP beads to reaction
2. Mix to homogeneous solution
3. Incubate at RT for 5 min
4. Briefly spin down
5. Place on magnet for 4 min
6. Discard supernatant carefully
7. Wash beads 2x with 200 μ L fresh 70% EtOH
8. Remove residual fluid with 20 μ L pipette
9. Air dry 3–5 min at RT
10. Elute DNA with volume* of 10 mM Tris-HCl pH 8.0
11. Place on magnet for 2 min
12. Remove purified product

* Refer to Product Insert for panel-specific parameters

† Work at room temperature for bead purifications

Appendix: Ligation cleanup bead purification†

1. Perform buffer exchange with 50 μ L of fresh Ligation cleanup buffer
2. Combine 50 μ L Ligation cleanup beads with entire volume Ligation Step 2 reaction
3. Mix solution by vortexing and incubate 5 min (2x)
4. Briefly spin down
5. Place tubes on magnet for 1 min
6. Discard supernatant carefully
7. Wash beads 2x with 200 μ L Ligation cleanup buffer
8. Wash beads with 200 μ L ultrapure H₂O
9. Resuspend beads in volume* 5 mM NaOH
10. Incubate reactions 10 min at 75 $^{\circ}\text{C}$
11. Cool reactions to 4 $^{\circ}\text{C}$
12. Briefly spin down
13. Place on magnet for 2 min
14. Remove purified product

For research use only. Not for use in diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.