

Primetime[™] qPCR Probes

Double- and single-quenched probes for use in your 5' nuclease assays



Dyes and quenchers for every experiment

PrimeTime qPCR Probes can be used in demanding applications such as multiplexing and digital PCR. PrimeTime qPCR Probes are available in a wide variety of dye-quencher combinations (**Table 1**) that are compatible with common qPCR instruments.

Consistent Results

All PrimeTime probes are HPLC purified, and then quality controlled via mass spectrometry, to deliver batch-to-batch consistency and minimize the need for troubleshooting.

Table 1. Commonly used fluorophores and quenchers.

Fluorophore		Emission wavelength (nm)	Quencher
	6-FAM*	520	ZEN/Iowa Black™ FQ
	SUN™∗	554	
	JOE™*	555	
	HEX*	555	
	MAX™∗	557	
	Cy® 3	564	Iowa Black RQ ⁺⁺
	ATTO [™] 550§	575	
	ROX	608	
	Texas Red® -X	617	
	ATTO 647N§	662	
	Cy 5 [¥]	668	
	Cy 5.5	706	Black Hole Quencher®-3¶

- * Probes with 6-FAM, SUN, JOE, MAX, or HEX fluorophores are also available as traditional, single-quenched probes with either Iowa Black FQ (license free) or Black Hole Quencher-1 (additional third-party licenses required for diagnostic use).
- Black Hole Quencher-2 (BHQ2) may also be used as a quencher (additional third-party licenses required for diagnostic use).
- Double-quenched probes available as a custom order.
- § ATTO-labeled probes available as a custom order.
- ¥ Cy 5 is also available as a singlequenched probe with BHQ2 (additional third-party licenses required for diagnostic use).
- Available as research use only.



RUO For Research Use Only. Not for use in diagnostic procedures.

idtdna.com



Enhance your assays with double-quenched probes

Enhance your assays and reduce background with ZEN or TAO Double-Quenched Probes. Our internal quenchers are 9 bases from the 5' fluorophore and work in combination with the 3' lowa Black Quencher (**Figure 1**).



Figure 1. Schematic of a PrimeTime qPCR 5' Nuclease Assay using a double-quenched probe that includes a dye, a ZEN or TAO internal quencher, and a 3' quencher.

With nearly 4 times lower background fluorescence (**Figure 2A**) and approximately 30% increased signal (**Figure 2B**), ZEN Double-Quenched Probes simply perform better. See data for TAO Double-Quenched Probes at **www.idtdna.com/qPCRprobes**.



Figure 2. Increase signal from ZEN Double-quenched probes. (A) ZEN probes provide greater dye quenching, producing lower background and, therefore, higher signal intensities than standard single-quenched probes (BHQ probes). (B) ZEN probes show earlier observed C_q values compared to BHQ single-quenched probes. Three replicate reactions with each probe type (40 bases long) were run with a gBlocks[®] Gene Fragment template (2 × 10⁵ copies) and PrimeTime[®] Gene Expression Master Mix (IDT) on the QuantStudio 7 qPCR instrument (ThermoFisher Scientific).

Quenching for long probes

Effective quenching for ZEN Double-Quenched Probes as long as 40 bases (**Figure 3**) means more effective designs, even for AT-rich targets.



Figure 3. Only ZEN Double-Quenched Probes maintain low background signal with increasing probe length. Probes of either 20 or 40 bases with 4 different quenchers run in 3 replicate reactions with each probe type run with a gBlock Gene Fragment template (2 × 10⁵ copies) and PrimeTime Gene Expression Master Mix (IDT) on the QuantStudio 7 qPCR instrument (ThermoFisher Scientific). Key: BHQ1= Black Hole Quencher-1 (Biosearch Technologies) and FQ = Iowa Black FQ (IDT).

Ordering Information

Visit www.idtdna.com/qPCRprobes to enter your sequence and choose modifications.

For more information, visit www.idtdna.com/qPCRprobes



For Research Use Only. Not for diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2025 Integrated DNA Technologies, Inc. All rights reserved. For specific trademark and licensing information, see www.idtdna.com/trademarks. Doc ID: RU022-1570_002 04/25