

# Alt-R™ Genome Editing Detection

> SEE WHAT MORE WE CAN DO FOR YOU AT [WWW.IDTDNA.COM](http://WWW.IDTDNA.COM).

NGS • CRISPR • Functional Genomics • Synthetic Biology • PCR & qPCR • Custom Oligos • Contract Manufacturing

## Table of contents

Alt-R Genome Editing Detection Kit—Consumables	3
Process CRISPR-Cas9-edited genomic DNA from cultured cells	4
Amplify genomic DNA and detect mutations	5
Form heteroduplexes for T7E1 digestion	6
Visualize T7E1 mismatch detection results	7
References	9

# Alt-R GENOME EDITING DETECTION KIT— CONSUMABLES

Kit component*	25 rxn kit	100 rxn kit	1000 rxn kit
T7 Endonuclease I (T7EI, 1 U/μL)	50 μL	200 μL	2000 μL
T7EI Reaction Buffer (10X)	50 μL	200 μL	2000 μL
Alt-R™ Control A (template/primer mix)	20 μL	20 μL	20 μL
Alt-R™ Control B (template/primer mix)	20 μL	20 μL	20 μL

\* PCR master mix and target-specific primers/probe sets are not included in the kits.

 **Note:** Store the Alt-R Genome Editing Kit at  $-20^{\circ}\text{C}$ .

Use the Alt-R Genome Editing Detection Kit to determine on-target genome editing and estimate editing efficiency using T7 endonuclease I (T7EI). In the T7EI assay, target genomic regions from CRISPR-modified cells are amplified by PCR. The PCR products are denatured and reannealed to allow heteroduplex formation between wild-type DNA and CRISPR–mutated DNA. Mutations are then detected using T7EI, which recognizes and cleaves mismatched DNA heteroduplexes. T7EI assay results are analyzed by visualizing cleavage products and full-length amplicons by gel or capillary electrophoresis.

For a quick, robust assessment of editing events in CRISPR experiments, we recommend using the T7EI assay, instead of alternative methods such as Surveyor® mismatch endonuclease assays or Sanger sequencing. The T7EI method is simple, fast, and provides clean electrophoresis results. The T7EI assay is also compatible with a broad range of PCR buffers and does not usually require purification of PCR products before digestion. Note that T7EI activity is sensitive to the DNA:enzyme ratio, as well as incubation temperature and time [1]. T7EI recognizes insertions and deletions of  $\geq 2$  bases that are generated by non-homologous end joining (NHEJ) activity in CRISPR experiments [2]. Because T7EI does not recognize 1 bp indels, T7EI underrepresents the total editing.

Design PCR primers that amplify your experimental target site and adjacent sequence. We recommend using a 600–1000 bp PCR amplicon with  $>100$  bp flanking the CRISPR cut site and with the CRISPR cut site off-center to allow fragment resolution by gel analysis or capillary electrophoresis. Test the PCR assay to determine the optimal annealing temperature to use with your samples and to verify that only the expected PCR product is synthesized. You can design the PCR assays using the PrimerQuest™ Tool at [www.idtdna.com/PrimerQuest](http://www.idtdna.com/PrimerQuest). Calculate the  $T_m$  of your primers at [www.idtdna.com/OligoAnalyzer](http://www.idtdna.com/OligoAnalyzer).

Required materials	Ordering information*
Alt-R™ Genome Editing Detection Kit	IDT (Cat # 1075931, 1075932, 1075933)
KAPA HiFi HotStart PCR Kit	Kapa Biosystems (Cat # KK2501)
Nuclease-Free Water	IDT (Cat # 11-04-02-01)
PCR primers specific for your target	IDT Custom PrimeTime™ PCR Primers ( <a href="http://www.idtdna.com/PrimerQuest">www.idtdna.com/PrimerQuest</a> )
Phosphate buffered saline (PBS)	General laboratory supplier
QuickExtract™ DNA Extraction Solution	Epicentre (Cat # QE09050)
(Option 1) Agarose	General laboratory supplier
(Option 2) Mutation Discovery Kit†	Advanced Analytical Technologies, Inc. (Cat # DNF-910-K1000T)

\* These are suggested sources for reagents used by the IDT R&D team when this protocol was written.

Individual components (e.g., the polymerase and buffer from the PCR kit) may be substituted with some optimization.

† For use on a Fragment Analyzer™ system (Advanced Analytical Technologies, Inc.).

Go to [www.idtdna.com](http://www.idtdna.com) for safety data sheets (SDSs) and certificates of analysis (COAs) for IDT products. Visit [www.idtdna.com/protocols](http://www.idtdna.com/protocols) to verify that you are using the most current version of this protocol.

# PROCESS CRISPR-Cas9-EDITED GENOMIC DNA FROM CULTURED CELLS

---

The volumes in **Steps 2** and **5** are optimized for confluent samples in 96-well plates. Some cell types and samples in larger wells will require larger volumes.

1. Wash CRISPR-Cas9-treated cells with 100  $\mu$ L of PBS.
2. Lyse cells by adding 50  $\mu$ L of QuickExtract™ DNA Extraction Solution.
3. Transfer cell lysate to appropriate PCR tubes or plate.
4. Vortex and heat in a thermal cycler at 65°C for 10 min, followed by 98°C for 5 min.
5. Add 100  $\mu$ L of Nuclease-Free Water to dilute the genomic DNA.
6. Vortex and spin down.

# AMPLIFY GENOMIC DNA AND DETECT MUTATIONS

1. Set up the PCR using template, primers, and components of the Alt-R Genome Editing Detection Kit and KAPA HiFi HotStart PCR Kit as follows:

Component	Amount: sample	Amount: Alt-R Control A	Amount: Alt-R Control B
Genomic DNA (Process CRISPR-Cas9–edited genomic DNA from cultured cells, <a href="#">step 6</a> )	4 $\mu$ L (~40 ng)	—	—
Forward primer	300 nM	—	—
Reverse primer	300 nM	—	—
Alt-R Control A (template/primer mix)	—	1 $\mu$ L	—
Alt-R Control B (template/primer mix)	—	—	2 $\mu$ L
KAPA HiFi Fidelity Buffer (5X)*	5 $\mu$ L (1X)	5 $\mu$ L (1X)	5 $\mu$ L (1X)
dNTPs	1.2 mM (0.3 mM each)	1.2 mM (0.3 mM each)	1.2 mM (0.3 mM each)
KAPA HiFi HotStart DNA Polymerase (1 U/ $\mu$ L)*	0.5 U	0.5 U	0.5 U
<b>Total volume</b>	<b>25 <math>\mu</math>L</b>	<b>25 <math>\mu</math>L</b>	<b>25 <math>\mu</math>L</b>

\* The reagents in the Alt-R Genome Editing Detection Kit are also compatible with other polymerases and buffers.

2. Run the PCR using the following cycling conditions:

Step	Temperature (°C)	Time (min:sec)	Cycles
Denature	95	5:00	1
Denature	98	0:20	
Anneal	Variable (primer specific)* 64–67 for Alt-R™ Control A and B	0:15	30
Extend	72	0:30	
Extend	72	2:00	1

\* Optimal annealing temperatures have been determined using KAPA HiFi HotStart DNA Polymerase. You may need to optimize for other polymerases. Annealing temperature for additional Alt-R controls: Alt-R CRISPR-Cas9 HPRT Primers, human and mouse = 67°C and Alt-R CRISPR-Cas9 HPRT Primers, rat = 64°C.

# FORM HETERODUPLEXES FOR T7EI DIGESTION

1. Combine the following in an appropriate PCR tube:

Component	Amount: Sample	Amount: Homoduplex control*	Amount: Heteroduplex control*
PCR (Amplify genomic DNA and Detect Mutations, <b>step 2</b> )	10 $\mu$ L experimental target or Alt-R HPRT control	10 $\mu$ L Control A	5 $\mu$ L Control A 5 $\mu$ L Control B
T7EI Reaction Buffer (10X)	2 $\mu$ L	2 $\mu$ L	2 $\mu$ L
Nuclease-Free Water	6 $\mu$ L	6 $\mu$ L	6 $\mu$ L
<b>Total volume</b>	<b>18 <math>\mu</math>L</b>	<b>18 <math>\mu</math>L</b>	<b>18 <math>\mu</math>L</b>

\* Alt-R Controls A and B are positive controls for the T7E1 assay that each contain template and PCR primers. Full length PCR products from Controls A and B contain the same sequences, except for a 6 bp deletion in Control B.

2. Heat and cool PCR products in a thermal cycler as follows:

Step	Temperature ( $^{\circ}$ C)	Time
Denature	95	10 min
Ramp 1	95–85	Ramp rate $-2^{\circ}$ C/sec
Ramp 2	85–25	Ramp rate $-0.3^{\circ}$ C/sec

3. Combine the following in a microcentrifuge tube for the T7EI digestion:

Component	Amount ( $\mu$ L)
PCR heteroduplexes (Form heteroduplexes for T7EI digestion, <b>step 2</b> )	18
T7 endonuclease I (1 U/ $\mu$ L)	2
<b>Total volume</b>	<b>20</b>


4. Incubate the T7EI reaction at 37 $^{\circ}$ C for 60 min.

# VISUALIZE T7EI MISMATCH DETECTION RESULTS

Visualize the digestion using one of the following methods:


- Use agarose gels.
- Dilute digestion with 150  $\mu$ L of 0.1X IDTE and run on a Fragment Analyzer™ system with the Mutation Discovery Kit. See **Figure 1** for representative results.

The expected amplicon and digested product sizes for Alt-R Controls A/B and the Alt-R CRISPR-Cas9 HPRT Positive Controls are listed in **Figures 1** and **2**, respectively.




**Figure 1. Sample data from Fragment Analyzer system—T7EI digestion containing PCR products from Alt-R Controls A and B.** PCR using template and primers in Controls A and B (Alt-R Genome Editing Detection Kit) were run using KAPA HiFi HotStart DNA Polymerase (Kapa Biosystems). PCR cycling conditions (annealing temperature, 64°C), denaturation/reannealing conditions, and T7EI digestion conditions were as described in this protocol. Digestion reactions were analyzed on a Fragment Analyzer system (Advanced Analytical Technologies, Inc.). Sample 1 contains Control A PCR products, while Sample 2 contains Control A and B PCR products. The trace (left) shows results from Sample 2. The gel image (right) shows results for Samples 1 and 2. RFU = relative fluorescence unit; LM = lower marker; UM = upper marker.

A. T7EI digestion—CRISPR-Cas9 editing (HPRT crRNA)



Cas9	Human	Mouse	Rat
Full-length (bp)	1083	647	930
Fragment 1 (bp)	827	512	599
Fragment 2 (bp)	256	135	331

B. T7EI digestion—CRISPR-Cas12a editing (HPRT crRNA)



Cas12a	Human	Mouse	Rat
Full-length (bp)	1083	647	930
Fragment 1 (bp)	578	400	657
Fragment 2 (bp)	505	147	273

**Figure 2. Sample data from T7EI digestion of Alt-R CRISPR HPRT control samples.** Genomic DNA from CRISPR-Cas9 (A) or Cas12a (B) edited human, mouse, and rat HPRT controls were PCR amplified, digested using T7 endonuclease I, and run on the Fragment Analyzer system. Reference standard bands at 5000 bp (upper marker) and 35 bp (lower marker) are used to align the gel and analyze the results. Estimated band sizes for the cut and uncut positive control amplicons are listed in the tables. Cell lines: HEK-293 (human), Hepa1-6 (mouse), and RG2 (rat).



# REFERENCES

---

1. Mean RJ, Pierides A, Deltas CC, et al. **Modification of the enzyme mismatch cleavage method using T7 endonuclease I and silver staining.** *Biotechniques*. 2004;36(5):758-760.
2. Vouillot L, Thelie A, Pollet N. **Comparison of T7EI and Surveyor mismatch cleavage assays to detect mutations triggered by engineered nucleases.** *G3 (Bethesda)*. 2015;5(3):407-415.

## Alt-R™ Genome Editing Detection Kit

Technical support: [www.idtdna.com/ContactUs](http://www.idtdna.com/ContactUs)

For more than 30 years, IDT's innovative tools and solutions for genomics applications have been driving advances that inspire scientists to dream big and achieve their next breakthroughs. IDT develops, manufactures, and markets nucleic acid products that support the life sciences industry in the areas of academic and commercial research, agriculture, medical diagnostics, and pharmaceutical development. We have a global reach with personalized customer service.

> SEE WHAT MORE WE CAN DO FOR YOU AT [WWW.IDTDNA.COM](http://WWW.IDTDNA.COM).

**For Research Use Only. Not for use in diagnostic procedures.** Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2022 Integrated DNA Technologies, Inc. All rights reserved. Trademarks contained herein are the property of Integrated DNA Technologies, Inc. or their respective owners. For specific trademark and licensing information, see [www.idtdna.com/trademarks](http://www.idtdna.com/trademarks).  
Doc ID: RUO22-1244\_001 11/22

