

Follow the instructions and incubations for each step. All mixing steps should be performed on ice. Pipette up and down 8 times or vortex to mix after re-suspending each lysphere and spin down prior to incubations and transfers. For incubations, use a lid heated $\geq 100^{\circ}\text{C}$, except where specified otherwise.

DNA Fragmentation	<p>Add 10-250 ng DNA to DNA Fragmentation reagent for a total volume of 50 μL.</p> <p>Incubate as indicated 🖐️</p>	<table border="1"> <thead> <tr> <th>Step</th> <th>Incubation Temperature</th> <th>Incubation Time</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>4°C</td> <td>1 min</td> </tr> <tr> <td>2</td> <td>37°C</td> <td>12 min</td> </tr> <tr> <td>3</td> <td>72°C</td> <td>20 min</td> </tr> <tr> <td>4</td> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table>	Step	Incubation Temperature	Incubation Time	1	4°C	1 min	2	37°C	12 min	3	72°C	20 min	4	4°C	Hold
	Step	Incubation Temperature	Incubation Time														
1	4°C	1 min															
2	37°C	12 min															
3	72°C	20 min															
4	4°C	Hold															
End Repair	<p>Transfer 50 μL DNA Fragmentation mixture to End Repair reagent.</p> <p>Incubate as indicated 🖐️</p>	<table border="1"> <thead> <tr> <th>Step</th> <th>Incubation Temperature</th> <th>Incubation Time</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>25°C</td> <td>30 min</td> </tr> <tr> <td>2</td> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table> <p>AMPure[™] XP clean-up (125 μL). Elute in 20 μL</p>	Step	Incubation Temperature	Incubation Time	1	25°C	30 min	2	4°C	Hold						
	Step	Incubation Temperature	Incubation Time														
1	25°C	30 min															
2	4°C	Hold															
Ligation Step 1	<p>Transfer 20 μL End Repair mixture to Ligation Step 1 reagent.</p> <p>Incubate as indicated 🖐️</p>	<table border="1"> <thead> <tr> <th>Step</th> <th>Incubation Temperature</th> <th>Incubation Time</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>37°C</td> <td>15 min</td> </tr> <tr> <td>2</td> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table> <p>AMPure[™] XP clean-up (50 μL). Elute in 42 μL</p>	Step	Incubation Temperature	Incubation Time	1	37°C	15 min	2	4°C	Hold						
	Step	Incubation Temperature	Incubation Time														
1	37°C	15 min															
2	4°C	Hold															
MBC Adapters	<p>Transfer 40 μL Ligation Step 1 mixture to the MBC adapters.</p>																
	<p>Transfer the entire volume of the MBC adapters to Ligation Step 2 reagent.</p> <p>Incubate as indicated with <u>unheated lid</u> 🖐️</p> <p>STOP Optional stopping point before purification. Store at -10°C to -30°C.</p>	<table border="1"> <thead> <tr> <th>Step</th> <th>Incubation Temperature</th> <th>Incubation Time</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>22°C</td> <td>5 min</td> </tr> <tr> <td>2</td> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table> <p>Ligation clean-up beads (50 μL). Elution volume varies*, 5 mM NaOH, 75°C 10 min</p>	Step	Incubation Temperature	Incubation Time	1	22°C	5 min	2	4°C	Hold						
Step	Incubation Temperature	Incubation Time															
1	22°C	5 min															
2	4°C	Hold															
Ligation Step 2																	

First PCR

Add GSP1 primers (volume varies*) to the First PCR HS/HGC reagent. Then add the entire volume of purified Adapter Ligation mixture and mix.

Incubate as indicated 🖐️



Optional stopping point after this step. Store at -10°C to -30°C.

Incubation Temperature	Incubation Time	# of cycles
95°C	3 min	1
95°C	30 sec	Varies*
Varies* (ramp rate 100%)		
72°C	3 min	1
4°C	Hold	1

AMPure™ XP clean-up (32 μL). Elution varies*

Second PCR

Add GSP2 primers (volume varies*) to the Second PCR HS/HGC reagent. Then add the entire volume of purified Second PCR reagent and mix.

Incubate as indicated 🖐️



Optional stopping point after this step. Store at -10°C to -30°C.

Incubation Temperature	Incubation Time	# of cycles
95°C	3 min	1
95°C	30 sec	Varies*
Varies* (ramp rate 100%)		
72°C	3 min	1
4°C	Hold	1

AMPure™ XP clean-up (32 μL). Elute in 20 μL

*Refer to Product Insert for panel-specific parameters

Proceed with Protocol: Quantify, Normalize, and Sequence Protocol for Illumina®

Appendix: AMPure® XP bead purification[†]

- Add AMPure XP beads to reaction
 - Mix to homogeneous solution
 - Incubate at RT for 5 min
 - Briefly spin down
 - Place on magnet for 4 min
 - Discard supernatant carefully
 - Wash beads 2x with 200 μL fresh 70% EtOH
 - Remove residual fluid with 20 μL pipette
 - Air dry 3-5 min at RT
 - Elute DNA with volume* of 10 mM Tris-HCl pH 8.0
 - Place on magnet for 2 min
- purified product

*Refer to Product Insert for panel-specific parameters

[†]Work at room temperature for bead purifications

Appendix: Ligation cleanup beads purification[†]

- Perform buffer exchange with 50 μL of fresh Ligation cleanup buffer
 - Combine 50 μL Ligation cleanup beads with 50 μL Ligation Step 2 reaction
 - Mix solution by vortexing and incubate 5 min (2x)
 - Briefly spin down
 - Place tubes on magnet for 1 min
 - Place on magnet for 4 min
 - Discard supernatant carefully
 - Wash beads 2x with 200 μL Ligation cleanup buffer
 - Wash beads with 200 μL ultrapure H₂O
 - Resuspend beads in volume* 5 mM NaOH
 - Incubate reactions 10 min at 75°C
 - Cool reactions to 4°C
 - Briefly spin down
 - Place on magnet for 2 min
- purified product