

# PrimeTime™ Gene Expression Master Mix

Product	Quantity*	Catalog #
PrimeTime Gene Expression Master Mix	1 x 1 mL	1055770
	1 x 5 mL	1055772
	5 x 5 mL	1055771
	20 x 25 mL	1072102

\* Separate tubes of reference dye are included with all orders, except for the 20 x 25 mL size (Reference dye can be ordered if needed. For more information visit, [www.idtdna.com/ContactUs](http://www.idtdna.com/ContactUs).)

Go to [www.idtdna.com/qPCRmastermix](http://www.idtdna.com/qPCRmastermix) for safety data sheets (SDSs) and certificates of analysis (COAs) for IDT products. Visit [www.idtdna.com/protocols](http://www.idtdna.com/protocols) to check that you are using the most current version of this protocol.

## Contents and storage conditions

PrimeTime Gene Expression Master Mix is a ready-to-use, 2X concentrated master mix that is designed for use in probe-based, real-time quantitative PCR. PrimeTime Gene Expression Master Mix contains an antibody-mediated, hot-start DNA polymerase, dNTPs, MgCl<sub>2</sub>, enhancers, and stabilizers.

In addition, a reference dye is provided as a separate component, making this master mix compatible for use on both reference dye-dependent and -independent instrument systems.

PrimeTime Gene Expression Master Mix is shipped at ambient temperature. Upon receipt, we recommend that you store the PrimeTime Gene Expression Master Mix at -15 to -30°C in a constant temperature (non-frost-free) freezer until the end of the month indicated in the expiration date, or at 2–8°C up to 1 month. We also recommend avoiding multiple freeze-thaw cycles and prolonged exposure of the reference dye to light.

## Protocol

### Prepare reagents

1. Thaw the PrimeTime Gene Expression Master Mix, reference dye, primer and probe solutions, and template DNA on ice.
2. Briefly vortex to mix reagents thoroughly.
3. Quickly centrifuge to collect solutions at the bottom of tubes.

### Add reference dye, as needed

If needed (see [Table 1](#)), add the appropriate amount of reference dye to the PrimeTime Gene Expression Master Mix ([Table 2](#)).

**Table 1. Reference dye concentration levels required by various PCR systems.**

PCR system*	Reference dye		
	High	Low	None
7900HT Fast and 7300 Real-Time PCR System (Thermo Fisher Scientific)	X		
StepOne™ and StepOnePlus™ Real-Time PCR System (Thermo Fisher Scientific)	X		
Mx3005P™ and Mx4000P™ qPCR System (Agilent)		X	
7500 Real-Time PCR System (Thermo Fisher Scientific)		X	
Viiia™ 7 Real-Time PCR System (Thermo Fisher Scientific)		X	
QuantStudio™ Flex Systems (Thermo Fisher Scientific)		X	
CFX, iQ™, and Opticon™ Real-Time PCR Detection Systems (Bio-Rad)			X
LightCycler® Real-Time PCR Systems (Roche)			X

\* For instruments not listed, please check with the manufacturer.

**Table 2. Amount of reference dye to add to PrimeTime® Gene Expression Master Mix.**

	Dye volume (µL)	
	High reference dye systems	Low reference dye systems
To a 1 mL stock vial of master mix, add:	40	4
To a 5 mL stock vial of master mix, add:	200	20
To a 25 mL stock vial of master mix, add:	1000	100

# Prepare the Reaction Mix

Make enough Reaction Mix for the number of reactions needed.

- Note:** Include controls (no reverse transcriptase, when needed; no template control; positive control) and 1–3 additional reactions to account for pipetting errors.

Combine and thoroughly mix all components (Table 3), except for the DNA template, which will be added separately in Add DNA template step.

- Note:** Add enough Nuclease-Free Water (Catalog # 11-04-02-01 or 11-05-01-14) so that the total reaction volume is 20  $\mu$ L or 10  $\mu$ L, depending on your chosen reaction volume.

Table 3. Reaction mixes.

Standard RT-qPCR reaction using purified RNA			
Component	Final concentration or amount	Volume per 20 $\mu$ L reaction**	Volume per 10 $\mu$ L reaction
PrimeTime Gene Expression Master Mix (2X)	1X	10 $\mu$ L	5 $\mu$ L
PrimeTime qPCR Assay* (20X)	1X	1 $\mu$ L	0.5 $\mu$ L
DNA template (Add in Add DNA template step)	3 pg to 100 ng	2–5 $\mu$ L	2–4.5 $\mu$ L
Nuclease-Free Water		Bring to 20 $\mu$ L	Bring to 10 $\mu$ L

  

qPCR using separate primers and probe			
Component	Final concentration or amount	Volume per 20 $\mu$ L reaction**	Volume per 10 $\mu$ L reaction
PrimeTime Gene Expression Master Mix (2X)	1X	10 $\mu$ L	5 $\mu$ L
Forward and reverse primers	250–1000 nM each	Variable	Variable
Probe(s)	150–250 nM each	Variable	Variable
DNA template (Add in Add DNA template step)	3 pg to 100 ng	2–5 $\mu$ L	2–4 $\mu$ L
Nuclease-Free Water		Bring to 20 $\mu$ L	Bring to 10 $\mu$ L

\* To order, visit [www.idtdna.com/PrimeTime](http://www.idtdna.com/PrimeTime).

\*\* The volumes provided are “per reaction”. Thus, determine final component volumes by multiplying each volume by the total number of reactions.

1. Dispense equal aliquots of the Reaction Mix into the wells of a qPCR plate that is compatible with your real-time PCR instrument.

## Add DNA template

1. Add 3 pg to 100 ng of DNA template or controls to the wells containing Reaction Mix (Table 3).
2. Seal with optically transparent film.
3. Gently agitate to make sure the mix is thoroughly combined, then briefly centrifuge to remove air bubbles and collect the reaction at the bottom of the wells.

## Set up cycling program

1. Program the appropriate PCR cycling protocol on your real-time PCR instrument (Table 4).

Table 4. Cycling protocol.

Step	Cycles	Temperature	Fast cycling (min:sec)	Standard cycling (min:sec)
Polymerase activation	1	95°C	3:00	3:00
Amplification:	35–45			
Denaturation		95°C	0:05	0:15
Annealing/Extension*		60°C	0:30	1:00
Hold, if needed	1	4°C	Up to 24 hr	Up to 24 hr

\* This is a general starting point. The annealing/extension temperature or time may need to be adjusted based on primer sequences. However, **do not change the polymerase activation conditions**.

## Run PCR

1. If the plate was stored before PCR, vortex, then briefly centrifuge.
2. Place the plate in the real-time PCR instrument and start the cycling program.

For additional information or assistance, go to [www.idtdna.com/ContactUs](http://www.idtdna.com/ContactUs).

**For Research Use Only. Not for use in diagnostic procedures.** Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2022 Integrated DNA Technologies, Inc. All rights reserved. Trademarks contained herein are the property of Integrated DNA Technologies, Inc. or their respective owners. For specific trademark and licensing information, see [www.idtdna.com/trademarks](http://www.idtdna.com/trademarks). Doc ID: RUO22-1161\_001 10/22