96-well plates with primers only

Primers (delivered amount: 5 nmol each) are shipped dried down and can be stored dried for up to 2 weeks at room temperature.

- 1. Centrifuge PrimeTime Primer Plates at 750 x g for 10 sec using a centrifuge with an appropriate plate adapter to ensure the product is at the bottom of the plate.
- 2. Add IDTE Buffer (Catalog # 11-05-01-13; 10 mM Tris, 0.1 mM EDTA, pH 8.0) to each well, pipetting the solution up and down the sides of each well to ensure maximum product recovery. Use the recommended volume indicated in the following table:

	Final desired stock concentration		
	40X	20X	10X
Well of PrimeTime™ Primer Plate	250 μL	500 μL	1000 μL



Note: PrimeTime qPCR Primers can be resuspended as 40X, 20X, or 10X stocks. Different concentrations may be preferable depending on the size of your PCR assay(s) and final reaction volume. The resuspended primers will yield a final 1X concentration of 500 nM primers.

- 3. Centrifuge the plate at 750 x g for 10 sec to collect the resuspended primers at the bottom of the wells.
- 4. Seal plate with a new adhesive plate cover and store resuspended primers at -20°C.

For Research Use Only. Not for use in diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2022 Integrated DNA Technologies, Inc. All rights reserved. Trademarks contained herein are the property of Integrated DNA Technologies, Inc. or their respective owners. For specific trademark and

© 2022 Integrated DNA Technologies, Inc. or their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks. Doc ID# RUO22-1058_001 08/22



