PrimeTime[™] qPCR Probe Assay

- 1. Centrifuge PrimeTime qPCR Probe Assay tubes at 750 x g for 10 seconds to ensure that primers and probes are at the bottom of the tube.
- 2. Resuspend the probe assay in IDTE buffer (Cat no. 11-01-02-05) to the recommended volume, depending on the size of the probe assay and the final reaction volume:

	Final stock concentration		
	40X	20X	10X
PrimeTime Mini qPCR Probe Assay (100 rxn)	Not recommended	100 µL	200 µL
PrimeTime Std qPCR Probe Assay (500 rxn)	250 µL	500 µL	1000 µL
PrimeTime XL qPCR Probe Assay (2500 rxn)	1250 µL	Not recommended	Not recommended

- **Note:** During resuspension, pipette mix the solution along the sides of the tube to ensure maximum product recovery. The number of reactions is based on a 20 μL reaction volume.
- 3. Centrifuge the resuspended the probe assay at 750 x g for 10 seconds.
 - Note: The resuspended assays will yield a final 1X concentration of 500 nM primers and 250 nM probe for probe assays ordered under default conditions. Standard or XL probe assays ordered with a custom primer:probe ratio will yield a different concentration of primers.

Storage: PrimeTime qPCR Probe Assays are shipped dried down and can be stored up to 2 weeks at room temperature, if shielded from light. Store resuspended assays in the dark at -20°C.

For Research Use Only. Not for use in diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations. © 2022 Integrated DNA Technologies, Inc. All rights reserved. Trademarks contained herein are the property of Integrated DNA Technologies, Inc. or their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks. Doc ID: RUO22-1225_001 12/22





WWW.IDTDNA.COM