Sequencing SARS-CoV-2

With the Midnight Amplicon Panel

This protocol uses:

xGen[™] DNA Library Prep Kit EZ

The method presented here is provided by IDT. This method can be used as a starting point for using SARS-CoV-2-Midnight-1200 Amplicon Panel, 500 rxn, and the xGen™ DNA Library Prep Kit EZ in similar experiments but may not be fully optimized for your application. IDT does not guarantee this method, as it has not been fully verified by our internal team. As a result, application scientists at IDT can only provide general guidance with limited troubleshooting support.

The products and protocol referenced herein are for research use only. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Any purchaser is solely responsible for all decisions regarding the use of these products and protocol and any associated regulatory or legal obligations.

Overview

COVID-19, caused by the coronavirus SARS-CoV-2, is an unprecedented global health challenge. The coronavirus has spread rapidly and causes a range of symptoms from mild upper respiratory effects to critical respiratory illness and even death. Understanding the biology of the virus is vital to combat it and reduce the clinical burden caused by this infection. Coronaviruses are RNA viruses and copy their genomic content directly, RNA to RNA. This is a relatively error-prone process, and over time mutations occur within the viral genome. Some of the mutations are neutral and randomly emerge and then disappear. Others provide some sort of advantage to the virus and become more common within the population. The emergence of new variants adds further challenges, leading to a need to track the virus.

Sequencing viral samples has been made more straightforward by amplification enrichment approaches that take total RNA from samples that include virus and human cells. Samples are reverse transcribed to form cDNA, and then specific amplification of the viral content is performed. To ensure complete viral genome coverage, two pools of primers are often used to separate overlapping amplicons. After amplification, the DNA is pooled and ligated to form an NGS-compatible library. Two sequencing platforms, Oxford Nanopore and Illumina®, are in common use for this process.

The Midnight Amplicon Panel was originally created by Drs Nikki Freed and Olin Silander of Massey University in New Zealand [1]. This panel consist of 29 amplicons of approximately 1200 bp in length. The reduced number of amplicons allows for more even sequencing across the genome than would be achievable with a higher number of amplicons. The panel was designed for use on the Oxford Nanopore sequencing platform. The length of the amplicons presents a problem for use on Illumina sequencers, which do not have read length options to cover the entire amplicon. This means that fragmentation of the amplicons is required before library construction. Here we present the use of the xGen DNA Library Prep Kit EZ as a simple workflow (Figure 1) to solve this challenge. The protocol includes a simple enzymatic fragmentation step which quickly creates lengths compatible with Illumina sequencing. Data indicate this workflow is highly efficient and rapid.

| 1 | cDNA synthesis | 30 minutes |
|---|---------------------------------|--------------------------------------|
| 2 | Amplicon generation | 2.5–3.5 hours |
| 3 | Post-amplification cleanup | 0.65X Ampure • 20 minutes |
| 4 | Fragmentation/End-repair/A-tail | Based on 350 bp fragment time in CoA |
| 5 | Ligation | 30 minutes |
| 6 | Post-ligation cleanup | 0.8X Ampure • 20 minutes |
| 7 | PCR amplification | 20–45 minutes |
| 8 | Post-PCR cleanup | 0.65X Ampure • 20 minutes |

Figure 1. Workflow of the demonstrated protocol.

Consumables and equipment

Consumables—IDT

| Item Catalog number | |
|---------------------------------------|---|
| xGen DNA Library Prep Kit EZ | www.idtdna.com/libraryprep 10009821 (96 rxn) |
| xGen Normalase™ UDI Primers Plate 1 | www.idtdna.com/NGS-adapters 10009796 |
| Nuclease-Free Water | 11-04-02-01 |
| Midnight Amplicon Panel, 500 rxn | website |
| 1 Liter IDTE, pH 8.0 (1X TE Solution) | 11-05-01-09 |

Consumables—Other suppliers

| Item | Supplier | Catalog number |
|---|--------------------------------|------------------|
| SuperScript™ IV | Thermo Fisher Scientific | 18091050 |
| Q5® High-Fidelity 2X Master Mix | New England BioLabs | M0492S |
| Absolute ethanol (200 proof) | General laboratory supplier | Varies |
| Agencourt® AMPure® XP–PCR purification beads, 5 mL | Beckman Coulter | A63880 |
| twin.tec® PCR plate 96 LoBind®, semi-skirted | Eppendorf | 0030129504 |
| Library Quantification Kit–Illumina/ Universal | Kapa Biosystems | KK4824 |
| MAXYMum Recovery® Microtubes, 1.7 mL | Corning | MCT-175-L-C |
| Microseal® 'B' PCR Plate Sealing Film | Bio-Rad | MSB1001 |
| Qubit™ dsDNA HS Assay Kit, or equivalent | Thermo Fisher Scientific | Q32851 or Q32854 |
| High Sensitivity DNA Kit, or equivalent | Agilent | 5067-4626 |
| High Sensitivity D1000 ScreenTape [™] , or equivalent | Agilent | 5067-5584 |

Equipment

| ltem | Supplier | Catalog number |
|--|-----------------------------|--------------------|
| 2200 TapeStation™ System/4200 TapeStation System, or equivalent | Agilent | G2965AA or G2991AA |
| Qubit 4 Fluorometer, or equivalent | Thermo Fisher Scientific | Q33226 |
| Microcentrifuge | General laboratory supplier | Varies |
| Plate centrifuge | General laboratory supplier | Varies |
| Thermal cycler | General laboratory supplier | Varies |
| qPCR system | General laboratory supplier | Varies |
| Magnet options (choose one): | | |
| Magnum™ EX Universal Magnet Plate | Alpaqua | A000380 |
| Magnetic Stand-96 | Thermo Fisher Scientific | Q33226 |

Template generation

The Midnight Amplicon Panel supports cDNA as input. Select a cDNA synthesis module [2] that supports 2 step RT-PCR, includes random primers, and has a processivity of >1 kb. This protocol has been tested with SARS-CoV-2 genomes isolated from nasopharyngeal swabs. Extracted RNA was used as input into the Superscript® IV First-Strand Synthesis System (Thermo Fisher Scientific Cat. No. 18091050), Lunascript™ RT Supermix kit (New England Biolabs Cat. No. E3010L) has also shown to be effective. The manufacturer's protocol was followed as written using the random primers and associated specifications. Optional host gDNA/RNA removal and RNAse H steps were not performed. The resulting cDNA was used directly as input into the Midnight Amplicon Panel.

cDNA synthesis

- 1. Thaw extracted RNA on ice.
- 2. Mix and briefly centrifuge all components prior to use.
- 3. Prepare the annealing Master Mix as follows on ice:

| Reagent | Volume per sample (μL) | |
|--------------------------|------------------------|--|
| 50 ng/µl random hexamers | 1 | |
| 10 mM dNTPs | 1 | |
| RNA | up to 11 | |
| DEPC-treated water | up to 11 | |
| Total volume | 13 | |

- 4. Seal, briefly vortex, and centrifuge samples.
- 5. Place samples on a thermal cycler and run the following program:

| Step | Temperature (°C) | Time |
|-----------|------------------|-----------|
| Annealing | 65 | 5 minutes |
| Hold | 4 | ∞ |

- 6. Place samples on ice for at least 1 min.
- 7. Make the RT Master Mix as follows:

| Reagent | Volume per sample (µL) | |
|------------------------|------------------------|--|
| 5xSSIV buffer | 4 | |
| 100 mM dNTPs | 1 | |
| Ribonuclease Inhibitor | 1 | |
| SSIV RT | 1 | |
| Total volume | 7 | |

- 8. Add 7 μ L RT Master Mix to each annealed reaction while on ice.
- 9. Seal, briefly vortex, and centrifuge samples.
- 10. Place samples on a thermal cycler and run the following program:

| Step | Temperature (°C) | Time |
|--------------|------------------|------------|
| Activation | 23 | 10 minutes |
| RT | 50 | 10 minutes |
| Inactivation | 80 | 10 minutes |
| Hold | 4 | ∞ |

Safe stopping point: Follow your established laboratory guidelines for proper storage of cDNA.

Amplicon generation

This portion of the protocol has been adapted for the Midnight Amplicon Panel from [2].

- 0
 - **Important:** Primer pools 1 & 2 are intended to be used separately and should not be combined. Each primer pool will need to be diluted before use.
 - 1. Prepare primer pools 1 & 2 by diluting the 100 μ M Midnight primer stocks to 10 μ M in IDTE.
 - 2. Mix and briefly centrifuge all components prior to use.
 - 3. Prepare Amplicon Master Mix as follows (one mix per primer pool):

| Reagent | Volume per sample (μL) | |
|-------------------------|------------------------|--|
| Q5 Hot Start Master Mix | 12.5 | |
| Midnight Primer Pool 1 | 1.1 | |
| Nuclease-Free Water | 8.9 | |
| Total volume | 22.5 | |

- 4. Mix thoroughly.
- 5. Add 2.5 μL of cDNA to the appropriate wells of a low-bind PCR plate (see **Note 1, Appendix A**). Each sample will need two wells (one well per primer pool).
- 6. Add 22.5 µL of Amplicon Master Mix to the appropriate wells.
- 7. Seal, briefly vortex, and centrifuge plate.
- 8. Place plate on a thermal cycler, and run the following program:

| Step | Number of cycles | Temperature (°C) | Time |
|-------------------------|---|------------------|------------|
| Heat Activation | 1 | 98 | 30 seconds |
| Denaturation | 25–35 (May need optimization; see Note 2, Appendix A) | 98 | 15 seconds |
| Annealing and Extension | | 65 | 5 minutes |
| Hold | 1 | 4 | ∞ |

- Safe stopping point: Follow your established laboratory guidelines for proper storage of cDNA.
- 9. Allow AMPure XP beads to equilibrate to room temperature. Vortex until homogenous before use.
- 10. Briefly centrifuge amplified cDNA plate.
- 11. Combine the entire volume of primer pool 1 and primer pool 2 PCR reactions per sample (total volume of 50μ L amplified cDNA) in a new low-bind PCR plate.
- 12. Add 32.5 µL (0.65X) of AMPure XP beads per sample.
- 13. Mix well by pipetting.
- 14. Incubate at room temperature for 5 minutes.
- 15. Transfer sample plate to a magnet and allow beads to form pellet (~2 minutes or until clear).

- 16. Remove and discard the supernatant without disturbing the bead pellet.
- 17. Wash the beads with 180 μ L of fresh 80% ethanol for 30 seconds, then remove the ethanol and discard. Do not disturb the pellet.
- 18. Repeat previous wash step for a total of two washes.
- 19. Using a new p20 pipette tip, remove any residual ethanol solution from the bottom of the well.
- 20. Allow beads to dry on the magnetic plate for 1–3 minutes.
- **Note:** Avoid over-drying beads, as this will reduce yield.
- 21. Remove sample plate from magnet, add 22 μL of Nuclease-Free Water, and resuspend by mixing.
- 22. Incubate at room temperature for 3 minutes.
- 23. Transfer sample plate back to the magnet, allow beads to collect for 5 minutes.
- 24. Carefully transfer 20 µL of supernatant to a new low-bind PCR plate, being careful to not disturb the pellet.
- **Safe stopping point:** Follow your established laboratory guidelines for proper storage of cDNA.

Note: To ensure an input of 1–250 ng into the library preparation, measure the concentration of the amplicon pool using a Qubit fluorometer (or equivalent).

Library preparation

Prepare reagents

- 1. Place all kit enzymes on ice (not in a cryocooler) for at least 20 minutes to allow enzymes to reach 4°C before pipetting.
- 2. Thaw other kit reagents, buffers, and primers on ice (4°C). Briefly vortex to mix well.
- 3. Spin all tubes in a microcentrifuge to collect contents before opening.
- 4. Prepare a fresh 80% ethanol solution using absolute ethanol (200 proof) and Nuclease-Free Water. At least 1.5 mL of 80% ethanol solution will be used per sample.

Perform enzymatic preparation

- Important: Keep the Enzymatic Prep Master Mix and the DNA samples on ice until they are loaded in the thermal cycler to safeguard against fragmentation; enzymes are active at room temperature and may fragment DNA to undesired sizes.
 - 1. Transfer the amplified cDNA sample (1–250 ng) to a new low-bind PCR plate, adjust the volume to a total of 19.5 µL (if necessary) using low EDTA TE Buffer, and place sample plate on ice.
 - **Note:** If not normalizing amplified cDNA concentration prior to library preparation, we recommend using 25 ng or more of input into the fragmentation.

| Reagent | Volume per sample (µL) | |
|----------------------|------------------------|--|
| Low EDTA TE 19.5 – x | | |
| Amplified cDNA | X | |
| Total volume | 19.5 | |

2. Set up the thermal cycler with the Enzymatic Prep program as described, with the heated lid set to 70°C:

| Step | Temperature (°C) | Time |
|------------------|------------------|------------|
| Initial hold | 4 | ∞ |
| F | 32 | See CoA |
| Enzymatic prep – | 65 | 30 minutes |
| Final hold | 4 | ∞ |



Note: See your **Certificate of Analysis (CoA)** for fragmentation time recommendations for the lot number you receive. Reaction times may need to be optimized for individual samples. Specifically, for sample inputs <25 ng, longer fragmentation may be required. To find the CoA, enter the lot number, which can be found on a label on the bottom of the kit, at **www.idtdna.com/COA**. We suggest using the 350 bp fragmentation time from the CoA.

The TapeStation trace (Figure 2) shows amplicon fragmentation using 16-, 18-, and 20-minute fragmentation times; the CoA recommendation for fragmentation to 200 bp using this lot of fragmentase was 16 minutes.




Figure 2. TapeStation trace showing effects of fragmentation times.

- 3. Begin Enzymatic Prep program to chill thermal cycler to 4°C.
- 4. Prepare the Enzymatic Prep Master Mix by adding the components in the order shown:

| Reagent | Volume per sample (μL) | | |
|--------------|------------------------|--|--|
| Buffer K1 | 3.0 | | |
| Reagent K2 | 1.5 | | |
| Enzyme K3 | 6.0 | | |
| Total volume | 10.5 | | |



Note: Keep reagents and master mix on ice as much as possible.

5. Gently vortex the Enzymatic Prep Master Mix for 5 seconds, then briefly centrifuge.



Important: Thoroughly mix the Enzymatic Prep Master Mix before and after adding to your DNA samples. Because this master mix is viscous, failure to mix thoroughly could result in incomplete fragmentation.

- 6. Add 10.5 μ L of premixed Enzymatic Prep Master Mix to each well containing DNA samples and low EDTA TE to reach a final volume of 30 μ L.
- 7. Thoroughly vortex sample tubes for 5 seconds to mix.
- 8. Briefly centrifuge the sample plate, then immediately place in the chilled thermal cycler, and advance the Enzymatic Prep program to the 32°C fragmentation step.
- 9. Continue the Enzymatic Prep program to completion.

Important: Fragmented samples can be kept at 4°C no more than 1 hour.

Perform Ligation

1. If you are using DNA inputs <25 ng, dilute the xGen Stubby Adapter in Duplex Buffer, as described:

| DNA input (ng) | Adapter dilution | Stock concentration (µM) |
|----------------|------------------|--------------------------|
| ≥25 | No dilution | 15 |
| 10 | 10-fold (1:10) | 1.5 |
| 1 | 20-fold (1:20) | 0.75 |

2. Prepare the Ligation Master Mix, adding the components in the order shown:

| Reagent | Volume per sample (μL) |
|------------------------------|------------------------|
| Low EDTA TE | 9 |
| Buffer W1 (room temperature) | 12 |
| Enzyme W3 | 4 |
| Reagent W5 (stubby adapter) | 5 |
| Total volume | 30 |



Important: Slowly pipette the viscous Buffer W1 to avoid bubbles and ensure accuracy. Keep reagents and master mix on ice except for the W1 ligation buffer.

- 3. Gently vortex Ligation Master Mix for 5 seconds.
- 4. When the Enzymatic Prep program is complete, add 30 μ L of premixed Ligation Master Mix to the plate containing your fragmented amplified cDNA samples (total volume is 60 μ L).
- 5. Thoroughly vortex to mix for 5 seconds.
- 6. Set up the thermal cycler with the Ligation program as described, with lid heating OFF:

| Step | Temperature (°C) | Duration (minutes) |
|----------|------------------|--------------------|
| Ligation | 20 | 20 |
| Hold | 4 | ∞ |

- 7. Run the samples in the thermal cycler with the Ligation program.
- 8. Immediately proceed to Clean up ligation reaction.

Clean up ligation reaction



Note: Make sure the AMPure XP beads are equilibrated to room temperature before you begin.

- 1. Vortex the beads until the solution is homogeneous.
- 2. Add 48 μ L of beads to each sample at room temperature (0.8X).
- 3. Vortex mix, then briefly centrifuge the sample plate.
- 4. Incubate the samples for 5 minutes at room temperature.
- 5. Place the sample plate on a magnet until the solution clears and a pellet has formed (~2 minutes).
- 6. Remove and discard the supernatant without disturbing the pellet (less than 5 µL may be left behind).
- 7. Add 180 μ L of freshly prepared 80% ethanol solution to the sample while it is still on the magnetic plate. Do not disturb the pellet.
- 8. Incubate for 30 seconds, then carefully remove the ethanol solution.
- 9. Repeat steps 7 and 8 for a second ethanol wash.
- 10. Using a new pipette tip, remove any residual ethanol solution from the bottom of the well.
- 11. Allow beads to dry on the magnetic plate for 1–3 minutes.

Note: Avoid over-drying beads, as this will reduce yield.

- 12. Add 20 μL of low EDTA TE Buffer to the sample plate.
- 13. Pipette mix until homogeneous.
- 14. Incubate at room temperature for 5 minutes.
- 15. Place the sample plate on a magnet for 2 minutes.
- 16. Transfer the supernatant containing your library to a new low-bind PCR plate, being careful to avoid any bead carryover.

Safe stopping point: Libraries can be stored overnight at -20°C.

Perform PCR and Cleanup

- 1. Add $5 \mu L$ of UDI Primer mix to each sample.
- 2. Add 25 μ L of the PCR Master Mix to the entire eluted library sample/UDI (25 μ L), then vortex mix (total volume is 50 μ L).
- 3. Briefly centrifuge the sample plate, then put on ice.
- 4. Set up the thermal cycler with the PCR program as described with a heated lid set to 105°C.

| Step | Cycles | Temperature (°C) | Duration (seconds) |
|-----------------|--------------------------|------------------|--------------------|
| Activate enzyme | 1 | 98 | 45 |
| Denature | | 98 | 15 |
| Anneal | Varies (see table below) | 60 | 30 |
| Extend | | 72 | 30 |
| Final extension | | 72 | 60 |
| Hold | | 4 | Hold |

| | Cycling recommendations for PCR-amplified workflows | | | |
|--------|--|------------------------------|--|-----|
| DNA | Minimum recommended cycles for >4 nM* with 350 bp insert | | Minimum recommended cycles for >500 ng* with 200 bp insert | |
| input | | xGen DNA Library Prep EZ UNI | xGen DNA Library Prep EZ | |
| 1 µg | 3** | 0–3 | 3** | 1–2 |
| 100 ng | 3** | 0–3 | 5 | 5 |
| 10 ng | 6–7 | 6–7 | 9 | 9 |
| 1 ng | 9–10 | 9–10 | 12 | 12 |
| 100 pg | 11–12 | 11–12 | 15 | 15 |

^{*} Additional PCR cycles are recommended to generate >12 nM for Normalase (Refer to the xGen DNA Library Prep EZ Kit protocol, Appendix A).

^{**}When indexing by PCR, a minimum of 3 cycles is required to attach adapter sequences, irrespective of whether a sufficient library amount is available following ligation.



Important: The number of cycles required to produce sufficient library for sequencing will depend on sample input quantity and quality. In the case of low-quality samples, the number of cycles required may vary based on the quality of the sample and amount of usable DNA present. Approximate guidelines for high-quality DNA are indicated in table above, but the exact number of cycles required must be empirically determined. Adjust cycle number as needed.

5. When the PCR program is complete, vortex the room-temperature AMPure XP beads until the solution is homogeneous.

- 6. Add 32.5 µL of AMPure XP beads (0.65X) to each sample well and mix by pipetting.
- 7. Incubate the sample plate for 5 minutes at room temperature.
- 8. Place the sample plate on a magnet until the solution clears and a pellet has formed (~2 minutes).
- 9. Remove and discard the supernatant without disturbing the pellet (less than 5 µL may be left behind).
- 10. Keeping the sample plate on the magnet, add 180 μL of freshly prepared 80% ethanol solution without disturbing the pellet.
- 11. Incubate for 30 seconds, then carefully remove the ethanol solution.
- 12. Repeat steps 10 and 11 for a second ethanol wash.
- 13. With a new pipette tip, remove any residual ethanol solution from the bottom of the well.
- 14. Allow beads to dry on the magnet for 1–3 minutes.

Note: Over-drying beads will reduce yield.

- 15. Add 22 μL of low EDTA TE Buffer to the sample wells and mix by pipetting until homogeneous.
- 16. Incubate samples at room temperature for 5 minutes.
- 17. Place the sample plate on a magnet for 2 minutes.
- 18. Transfer the supernatant (20 μ L) containing the final library to a new low-bind PCR plate. Be careful to avoid any bead carryover.

Safe stopping point: Libraries can be stored overnight at -20°C.




Figure 3. An example TapeStation trace of final library product for reference.

Sequencing

Libraries are now ready for sequencing on an Illumina platform via the manufacturer's methods. Details are available at here.

Appendix A

Note 1

In an independent experiment, we tried adding 6 μ L of cDNA into the amplification reaction instead of 2.5 μ L to investigate whether we would see an increase in amplified material; we found that the increase in cDNA volume in the amplification reaction did not increase final yield of the pools. Therefore, only amplicon pools generated using 2.5 μ L of cDNA were carried through the entire protocol.

Note 2

As stated in the **Amplicon Generation** section, cycling parameters may need optimization. Here, we adjusted cycle number based on Ct value. Examples are shown in **Data Table 1**.

Data Table 1. Example data of RNA sample Ct values and corresponding post-amplification yields.

| C _t value based on N1 assay from 2019-nCoV RUO Kit | Number of PCR cycles for amplicon generation | Yield after clean-up (ng/μL) |
|---|--|---------------------------------|
| 29.3 | 30 | ~1.8 |
| 31.4 | 35 | ~9.8 |
| 33.2 | 35 | ~5.8 |

References

- Freed NE, Vlková M, Faisal MB, Silander OK. Rapid and inexpensive whole-genome sequencing of SARS-CoV-2 using 1200 bp tiled amplicons and Oxford Nanopore Rapid Barcoding. Biol Methods Protoc. 2020;5(1):bpaa014.
 Published 2020 Jul 18. doi:10.1093/biomethods/bpaa014
- 2. Freed N and Silander O. nCoV-2019 sequencing protocol (RAPID barcoding, 1200bp amplicon) v4 (protocols. io.bh7hj9j6). protocolsio. Published 2020 Jul 02.

Sequencing SARS-CoV-2 with the Midnight Amplicon Panel

Technical support: applicationsupport@idtdna.com

For more than 30 years, IDT's innovative tools and solutions for genomics applications have been driving advances that inspire scientists to dream big and achieve their next breakthroughs. IDT develops, manufactures, and markets nucleic acid products that support the life sciences industry in the areas of academic and commercial research, agriculture, medical diagnostics, and pharmaceutical development. We have a global reach with personalized customer service.

> SEE WHAT MORE WE CAN DO FOR YOU AT WWW.IDTDNA.COM.

For Research Use Only. Not for diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2022 Integrated DNA Technologies, Inc. All rights reserved. Trademarks contained herein are the property of Integrated DNA Technologies, Inc. or their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks. Doc ID: Doc ID: RUO22-0867_001 06/22