VARIANTPlex™-HT Dx Liquid Reagents

Intended Use

For *in vitro* diagnostic use.

For professional use only.


The VARIANTPlex™-HT Dx is a set of reagents and consumables used to prepare sample libraries, using manual or automated methods, from Nucleic Acid (DNA or TNA) extracted from fresh-frozen or formalin-fixed, paraffin-embedded tissue or cells. User-supplied oligonucleotide panels are required for the preparation of libraries targeting specific genomic regions of interest. The generated sample libraries incorporate P5 and P7 sequences for downstream use.

Summary and Principles of Procedure

Library preparation utilizing Anchored Multiplex PCR (AMP) is a rapid and scalable method to generate target-enriched libraries for NGS. AMP technology can be used for applications in targeted DNA sequencing to generate a sequencing library. Designed for low nucleic acid input, this process delivers robust performance.

AMP utilizes unidirectional gene-specific primers (GSPs) that enrich for both known and unknown mutations. Adapters that contain both molecular barcodes and sample indices permit quantitative multiplex data analysis, read deduplication, and accurate mutation calling.

Workflow Overview

* First and Second PCR times vary based on the specific GSP panel as well as individual lab cycling conditions. See Product Insert for more information.

Recommendations

- Read through the entire protocol before starting your library preparation.
- Take note of safe stopping points throughout the protocol where samples can be safely frozen (-30° C to -10° C) to plan your workflow.
- Wipe down workstation and pipettes with nuclease and nucleic acid cleaning products (e.g., RNase AWAY, Thermo Fisher Scientific).
- If utilizing 96-well plates, make sure they are nuclease free, wells are of sufficient volume and plate seals are sufficient to prevent moisture loss.

Warnings and Precautions

- Use good laboratory practices to prevent contamination of samples by PCR products.
- Use nuclease-free PCR tubes, microcentrifuge tubes, and aerosol-barrier pipette tips.

- Verify that the thermal cycler used for library preparation is in good working order and currently calibrated according to manufacturer specifications. When using a thermal cycler with a 100% ramp rate >6°C/sec, set the ramp rate to go no higher than 6°C/sec for First PCR and Second PCR.
- Ensure Reaction cleanup beads (used repeatedly throughout the workflow) are equilibrated to room temperature (20°C to 25°C) and fully resuspended prior to use.
- Ligation Step 2 Components are viscous and must be fully mixed prior to aliquoting.
- Do not use reagents beyond labeled expiry date.
- Immediately report any serious incidents related to this product to IDT and the Competent Authority or other local regulatory authority in which the user and the patient reside.

Limitations

- For *in vitro* diagnostic use.
- For use with input masses of 10 to 200ng
- The enzymes provided are sensitive to frequent temperature changes and potential freeze-thaw events.
- Enzymes will not freeze when stored at recommended storage temperature range but may freeze during shipment.
- For best results, it is recommended that you store the enzymes provided in a -20°C benchtop cooler box and use this to transport enzymes between workspaces and freezers.
- Enzyme components will remain stable for up to 5 freeze-thaw cycles without effecting the functionality.
- Buffers will freeze at recommended storage temperature and should be thawed on ice prior to use. For best results, Buffers should be aliquoted to appropriate batch sizes for single use with adequate overage, such that aliquots are not re-frozen after thawing. Input nucleic acid (DNA) in EDTA-free buffer (pH 7-8) or ultrapure water is the optimal starting template for AMP library preparation. Do NOT use EDTA-containing buffers.
- Contact technical support (archer-tech@idtdna.com) for commercially available extraction kit recommendations.
- Use the maximum allowable input mass (200 ng) whenever possible. Higher input quantities enable more sensitive variant detection.

Reagents Preparation

- The volumes shown below need to be increased if you are preparing more than 24 libraries.
- Make at least **10 mL** fresh 10 mM Tris-HCl, pH 8.0 from 1M Tris-HCl, pH 8.0 and ultrapure water.
 - Mix **100 µL** 1M Tris-HCl, pH 8.0 with **9,900 µL** ultrapure water.
 - 10 mM Tris-HCl, pH 8.0 can be used for up to one week after mixing.
- Make at least **50 mL** fresh 70% ethanol from 100% ethanol and ultrapure water.
 - Add **35 mL** 100% ethanol to **15 mL** ultrapure water.
 - Tightly close the cap to minimize evaporation when not in use.
 - 70% ethanol is appropriate for use for up to one week after mixing.
- Make at least **1 mL** fresh 5 mM NaOH working stock from concentrated NaOH and ultrapure water.
 - If working from 1 M NaOH, add **5 µL** of 1 M NaOH to **995 µL** of ultrapure water to yield 5 mM final NaOH.
 - If working from 5 M, add **10 µL** of 5 M NaOH to **990 µL** of ultrapure water to yield 50 mM NaOH. Mix well and briefly spin down. Take **100 µL** of 50 mM NaOH and combine with **900 µL** of ultrapure water to yield 5 mM NaOH. Mix well and briefly spin down.

Product Components

The VARIANT[™]Pl[™]ex[™]-HT Dx Kit components are listed as follows:

VARIANT[™]Pl[™]ex[™]-HT Dx Liquid Reagents – 96 reactions, REF 10029097

VARIANT[™]Pl[™]ex[™]-HT Dx Liquid Reagents – 24 reactions, REF 10029100

Archer[™] Cleanup Reagents - 96 reactions, REF 10029134

Archer[™] Cleanup Reagents - 24 reactions, REF 10029142

Archer[™] Liquid Adapter Dx Kit Set A – 96 reactions, REF 10029148

Archer[™] Liquid Adapter Dx Kit Set B – 96 reactions, REF 10029151

Materials Provided

VARIANTPlex™-HT Dx Liquid Reagents (96 reaction, 10029097; 24 reaction, 10029100)

Product	Contents	REF
DNA Fragmentation Buffer	77 or 308μL of reaction buffer for enzymatic DNA shearing (-30°C to -10°C)	10029101, 10029098
DNA Fragmentation Enzyme	116 or 416μL of enzyme for cleaving DNA (-30°C to -10°C)	10029102, 10029099
Complete End Repair Buffer	213 or 851μL of reaction buffer for end polishing (-30°C to -10°C)	10029070, 10029088
Complete End Repair Enzyme A	24 or 96μL aqueous solution containing glycerol and enzyme (-30°C to -10°C)	10029068, 10029086
Complete End Repair Enzyme B	3.4 or 13.5μL aqueous solution containing glycerol and enzyme (-30°C to -10°C)	10029069, 10029087
Ligation Step 1 Buffer	89 or 356μL of reaction buffer for Ligation Step 1 (-30°C to -10°C)	10029072, 10029090
Ligation Step 1 Enzyme	7.2 or 29μL aqueous solution containing glycerol and enzyme (-30°C to -10°C)	10029071, 10029089
Ligation Step 2 Buffer	591 or 1190μL of reaction buffer for Ligation Step 2 (-30°C to -10°C)	10029074, 10029092
Ligation Step 2 Enzyme	9.6 or 38.5μL aqueous solution containing glycerol and enzyme (-30°C to -10°C)	10029073, 10029091
First PCR Buffer 2.0	269 or 1076μL of reaction buffer for first PCR (-30°C to -10°C)	10029076, 10029094
Second PCR Buffer 2.0	269 or 1076μL of reaction buffer for second PCR (-30°C to -10°C)	10029077, 10029095
PCR Enzyme	38.4 or 154μL aqueous solution containing glycerol and enzyme (-30°C to -10°C)	10029075, 10029093

Archer™ Cleanup Reagents (24 reaction, 10029142; 96 reaction, 10029134)

Description	Contents and Storage Temperature	REF
Ligation Cleanup Beads	1200μL tube(s) of buffered magnetic bead slurry for capture of ligated products (2°C to 8°C)	10029135
Ligation Cleanup Buffer	10.8 or 43.2mL of buffered solution for washing ligation cleanup binding reactions (2°C to 8°C)	10029144, 10029136
Reaction Cleanup Beads	11.8 or 35.2mL of buffered magnetic bead slurry for immobilization of nucleic acids (2°C to 8°C)	10029145, 10029137

Archer™ Liquid Adapter Dx Kit (Set A 10029148; Set B 10029151)

Description	Contents and Storage Temperature	REF
Liquid P5 Adapter Plate A (Index 1-96)	96-well plate containing 2μL of MBC Adapter oligo (-30°C to -10°C)	10029146
Liquid P7 Index Plate A (Index 1-96)	96-well plate containing 4μL second PCR indexing primer (-30°C to -10°C)	10029147
Liquid P5 Adapter Plate B (Index 97-192)	96-well plate containing 2μL of MBC Adapter oligo (-30°C to -10°C)	10029149
Liquid P7 Index Plate B (Index 97-192)	96-well plate containing 4μL second PCR indexing primer (-30°C to -10°C)	10029150

Materials Required but Not Supplied

- 1 M Tris-HCl, pH 8.0 (molecular biology grade)
- Ultrapure water (molecular biology grade)
- 200 mM Tris-HCl, pH 7.0 (for sequencing)
- 100% ethanol (ACS grade) Various
- Concentrated NaOH solution (ACS grade)
- RNase AWAY™ (Thermo Fisher Scientific-REF:7003)
- KAPA Universal Library Quantification Kit (KAPA Biosystems-REF:KK4824)
- Standard PCR thermal cycler
- Real-Time PCR thermal cycler
- qPCR tubes
- 0.2 mL PCR tubes or 96-well plates
- DynaMag™-96 Side Magnet (Thermo Fisher Scientific-REF:12331D)
- Microcentrifuge
- Plate centrifuge
- Pipettes (P10, P20, P200 and P1000) (Pipetman or equivalent)
- Sterile, nuclease-free aerosol barrier pipette tips
- Vortex mixer
- PCR tube cooling block
- Gloves
- Qubit® Fluorometer (3.0 or higher) (Thermo Fisher Scientific-REF:Q33216)
- Qubit dsDNA HS Assay Kit (Thermo Fisher Scientific-REF:Q32854)
- -20°C 1.5mL Tube Benchtop Cooler Box Various
- Gene-Specific Primers (compatible with Anchored Multiplex PCR)

Summary of Performance Characteristics

Library yield of 4nM or greater is expected when input masses are from 10-200ng. The table below summarizes the performance results for library yield of VARIANTPlex™-HT Dx. All results pass the acceptance criteria.

Primer Panel	Input Description (Manufacturer and part number)	Input mass (ng)	Replicates	Result (Pass criteria >90% OPA)
Gene Specific Primer Oligonucleotide Pool 3	Mimix Myeloid Cancer Panel (Horizon Discovery HD829)	10	3	Pass
		50	3	
		200	3	
	Water (Sigma-Aldrich W4502-50ML)	0	5	Pass
Gene Specific Primer Oligonucleotide Pool 4	Mimix™ Quantitative Multiplex, FFPE Reference Standard IVD (Horizon Discovery HD200-IVD)	10	3	Pass
		50	3	
		200	3	
	Water (Sigma-Aldrich W4502-50ML)	0	5	Pass

Procedure

Before beginning, review “Working with Liquid Reagents” section for additional information regarding best practices for enzymes and buffers. For ease begin by programming thermal cycler to include the following:

Incubation Program	Step	Temperature (°C)	Time (min)
DNA Fragmentation	1	4	1
	2	30	4
	3	65	30
	4	4	Hold
Second Strand cDNA Synthesis	1	16	60
	2	75	20
	3	4	Hold
Complete End Repair (Heated lid off)	1	25	30
	2	4	Hold
Ligation Step 1	1	37	15
	2	4	Hold
Ligation Step 2 (Heated lid off)	1	22	5
	2	4	Hold

Ligation Elution	1	75	10
	2	4	Hold

PCR Cycling Programs	Step	Temperature (°C)	Time	Cycles
First PCR	1	95	3 min	1
	2	95	30 sec	Varies—See panel specific product insert)
	3	Varies—See panel specific product insert)	10 sec	
	4		Varies—See panel specific product insert) 100% ramp rate*	
	5	72	3 min	1
	6	4	Hold	1
Second PCR	1	95	3 min	1
	2	95	30 sec	Varies—See panel specific product insert)
	3	Varies—See panel specific product insert)	10 sec	
	4		Varies—See panel specific product insert) 100% ramp rate*	
	5	72	3 min	1
	6	4	Hold	1

*See ramp rate guidelines in the "Warnings and Precautions" section.

Step 1: DNA Fragmentation

1. Thaw the **DNA Fragmentation Buffer** on ice. Briefly vortex Buffer, then spin down and return to ice. **Do not vortex the enzyme**, instead gently invert the tube 3-5 times and flick to mix, then spin down.
2. Adjust purified DNA samples to a final volume of **32 µL** and transfer to a new 0.2 mL 8-strip tube or plate wells. Keep samples on ice.

Component	Reaction Mix
Ultrapure water	32 - Xµl
Purified nucleic acid	Xµl
Total volume	32µl

3. On ice, make a **DNA Fragmentation Master Mix** working solution based on the desired number of reactions (including but not exceeding 20% overage) using the table below.

DNA Fragmentation Master Mix	1X (µL)	___X (µL)
DNA Fragmentation Buffer (10029101 or 10029098)	3.2	
DNA Fragmentation Enzyme (10029102 or 10029099) Do NOT vortex, invert 3-5x and flick tube before adding	4.8	
Total volume	8	

- a. Mix by briefly vortexing and spin down.
 - b. Keep master mix on ice.
4. To each **32 µL** input sample, add **8 µL** DNA Fragmentation Master Mix.
 - a. Mix by briefly vortexing and spin down.
 - b. Return tubes to ice.
 5. Start the following thermal cycler program, and only transfer reactions to the block once temperature reaches 4°C. Pause the program if necessary.
 - a. Use a heated lid (≥100°C).

DNA Fragmentation incubation conditions

Step	Temperature (°C)	Time (minutes)
1	4	1
2	30	4
3	65	30
4	4	Hold

- b. After the program has reached 4°C, briefly spin down reactions and place on ice.

Step 2: Complete End Repair

- Thaw the **Complete End Repair Buffer** on ice. Briefly vortex Buffer and either vortex or pipette mix enzymes, then spin down.
- On ice, make a **Complete End Repair Master Mix** working solution based on the desired number of reactions (including but not exceeding 20% overage) using the table below.

Complete End Repair Master Mix	1X (μL)	___X (μL)
Complete End Repair Buffer (10029070 or 10029088)	8.86	
Complete End Repair Enzyme A (10029068 or 10029086)	1	
Complete End Repair Enzyme B (10029069 or 10029087)	0.14	
Total Volume	10	

- Mix by briefly vortexing and spin down.
 - Keep master mix on ice.
- To each **40 μL** DNA sample from Step 1: DNA Fragmentation, add **10 μL** of Complete End Repair Master Mix.
 - Mix by briefly vortexing and spin down.
 - Return tubes to ice.
 - Transfer reactions to a preheated thermal cycler and initiate an incubation using the following program and guidelines:
 - Heated lid off. If the lid is still hot, consider using a different thermal cycler or leaving the lid open.

Complete End Repair incubation conditions

Step	Temperature (°C)	Time (minutes)
1	25	30
2	4	Hold

- Place samples in the thermal cycler and start the program.
- When the run has completed, briefly spin down reactions and place on ice.

Reaction Cleanup after Complete End Repair

Refer to **Warnings and Precautions** section of this Instructions for Use for guidance on working with Reaction Cleanup beads.

- Completely resuspend Reaction Cleanup beads by vortexing.
- Remove tubes from ice and add **2.5X** volume (**125 μL**) of Reaction Cleanup beads to each Complete End Repair reaction.
- Vortex well or pipette 10 times to mix and visually inspect the color of the sample to ensure a homogenous mixture.
- Incubate for **5 minutes** at room temperature (20°C to 25°C).
- Briefly spin down tubes.
- Place tubes on the magnet for **4 minutes or until beads are fully pelleted** against the tube wall.
- Without disturbing the bead pellet, use a pipette to remove and discard the supernatant. If the pellet becomes dislodged from the magnet and a portion is drawn into the pipette tip, return contents to tube and repeat magnet incubation step.
- Wash beads **two times** with 70% ethanol while still on the magnet. For each wash:
 - Add **180 μL** 70% ethanol.
 - Incubate for **30 seconds** at room temperature (20°C to 25°C).
 - Carefully remove ethanol and discard.
- After the final wash, use a pipette (≤20 μL capacity) to completely remove visible supernatant residue and allow tubes to dry for **3-5 minutes** at room temperature with open lids. **Take care not to over-dry beads** as this will significantly decrease overall recovery (yield) of nucleic acid.
- Elute DNA by resuspending beads in **16 μL** 10mM Tris-HCl, pH 8.0.

***Note that for automated workflows it is acceptable to add 1 μL 10mM Tris-HCl, pH 8.0 to the bead resuspension volume shown above. This extra 1 μL should be left behind when removing the purified DNA in upcoming steps.**
- Place tubes back on the magnet for **2 minutes**.

Step 3: Ligation Step 1

1. Thaw the **Ligation Step 1 Buffer** on ice. Briefly vortex Buffer and either vortex or pipette mix enzymes, then spin down.
2. On ice, make a **Ligation Step 1 Master Mix** working solution based on the desired number of reactions (including but not exceeding 50% overage) using the table below.

Ligation Step 1 Master Mix	1X (μL)	___X (μL)
Ligation Step 1 Buffer (10029090 or 10029072)	3.7	
Ligation Step 1 Enzyme (10029089 or 10029071)	0.3	
Total Volume	4	

- a. Mix by briefly vortexing and spin down.
- b. Keep master mix on ice.
3. Transfer **16 μL** of eluted DNA from Step 2: Complete End Repair into new 0.2mL 8-strip tubes or plate wells. It is acceptable for a small amount of Reaction Cleanup beads to be transferred.
4. Add **4 μL** of Ligation Step 1 Master Mix to each sample.
- a. Mix by briefly vortexing and spin down.
- b. Return tubes to ice.
5. Transfer reactions to a preheated thermal cycler and initiate an incubation using the following program and guidelines:
 - a. Use a heated lid ($\geq 100^{\circ}\text{C}$).

Ligation Step 1 incubation conditions

Step	Temperature ($^{\circ}\text{C}$)	Time (minutes)
1	37	15
2	4	Hold

- b. Place samples in the thermal cycler and start the program.
- c. After the program has reached 4°C , remove tubes from the temperature block, briefly spin down reactions and place on ice.

Reaction Cleanup after Ligation Step 1

Refer to **Warnings and Precautions** section of this Instructions for Use for guidance on working with Reaction Cleanup beads.

1. Completely resuspend Reaction Cleanup beads by vortexing.
2. Remove tubes from ice and add **2.5X** volume (**50 μL**) of Reaction Cleanup beads to each Ligation Step 1 reaction.
3. Vortex well or pipette 10 times to mix and visually inspect the color of the sample to ensure even mixing.
4. Incubate for **5 minutes** at room temperature (20°C to 25°C).
5. Briefly spin down tubes.
6. Place tubes on the magnet for **4 minutes or until beads are fully pelleted** against the tube wall.
7. Without disturbing the bead pellet, use a pipette to remove and discard the supernatant. If the pellet becomes dislodged from the magnet and a portion is drawn into the pipette tip, return contents to tube and repeat magnet incubation step.
8. Wash beads **two times** with 70% ethanol while still on the magnet. For each wash:
 - a. Add **180 μL** 70% ethanol.
 - b. Incubate for **30 seconds** at room temperature (20°C to 25°C).
 - c. Carefully remove ethanol and discard.
9. After the final wash, use a pipette (≤ 20 μL capacity) to completely remove visible supernatant residue and allow tubes to dry for **3-5 minutes** at room temperature with open lids. Take care not to over-dry beads as this will significantly decrease overall recovery (yield) of nucleic acid.
10. Elute DNA by resuspending beads in **25 μL** 10 mM Tris-HCl, pH 8.0.
11. Place tubes back on the magnet for **2 minutes**.

Step 4: MBC Adapter Incorporation

1. Thaw **Liquid P5 Adapter Plate(s)** (10029146 or 10029149). Spin down before unsealing plate. It is acceptable to pierce the foil seal with a pipette tip to access the liquid adapters.

Important As this step incorporates the **P5 index tag** for sample-level tracking, be sure to record which **MBC adapter** is being used for each sample.
2. Transfer **23 μL** of eluted DNA from the cleanup after Step 3: Ligation Step 1 above to new 0.2mL 8-strip tubes or plate wells. **Avoid pipetting Reaction Cleanup beads into this reaction.**

3. Add **2 μ L** of respective unique **Liquid P5 MBC Adapter** for a total of **25 μ L**. Re-seal adapter plate using a new seal after all P5 MBC adapters have been added to respective samples.
 - a. Mix by vortexing and spin down.
 - b. Return tubes to ice.
4. Immediately proceed to Step 5: Ligation Step 2.

Step 5: Ligation Step 2

1. Thaw the **Ligation Step 2 Buffer** on ice. Vortex Buffer well and spin down. Either vortex or pipette mix enzymes then spin down.
2. On ice, make a **Ligation Step 2 Master Mix** working solution based on the desired number of reactions (including but not exceeding 20% overage) using the table below.

! Due to the viscosity of the Ligation Step 2 Buffer components, it is essential to slowly pipette mix the stock solution at least 10 times then vortex for 10 seconds BEFORE aliquoting the required volume for the working solution and again AFTER the addition of the Ligation Step 2 Enzyme.

Ligation Step 2 Master Mix	1X (μ L)	___X (μ L)
Ligation Step 2 Buffer (10029092 or 10029074)	24.6	
Ligation Step 2 Enzyme (10029091 or 10029073)	0.4	
<i>Total volume</i>	25	

- a. Mix by briefly vortexing and spin down.
 - b. Keep master mix on ice.
3. Transfer **25 μ L** of Ligation Step 2 Master Mix into each tube or well containing the **25 μ L** of each sample (DNA + Liquid MBC Adapter) from Step 4: MBC Adapter Incorporation above.
 - a. Mix well and spin down. **Due to the viscosity of these components, carefully pipette mixing 10 times in addition to vortexing is required to ensure complete mixing.**
 - b. Return tubes to ice.
4. Transfer reactions to a preheated thermal cycler and initiate an incubation using the following program and guidelines:
 - a. Heated lid off

Ligation Step 2 incubation conditions

Step	Temperature ($^{\circ}$ C)	Time (minutes)
1	22	5
2	4	Hold

- b. After the program has reached 4° C, briefly spin down reactions and place on ice.

— Safe stopping point: It is okay to stop and store the reactions at -30° C to -10° C.

Reaction Cleanup after Ligation Step 2

! Caution: This step uses Ligation Cleanup Beads (10029135 or 10029135) and Ligation Cleanup Buffer (10029144 or 10029136) instead of Reaction Cleanup beads and 70% ethanol. Make sure to use Ligation Cleanup Beads and Ligation Cleanup Buffer that come with the kit for this step.

Prepare Ligation Cleanup Beads:

1. Completely resuspend **Ligation Cleanup Beads** by vortexing.
2. For each reaction, pipette **50 μ L** of Ligation Cleanup Beads into new 0.2 mL 8-strip tubes.
3. Place tube(s) on the magnet for **1 minute or until the beads are pelleted**.
4. Without disturbing the bead pellet, use a pipette to remove and discard the supernatant. If the pellet becomes dislodged from the magnet and a portion is drawn into the pipette tip, return contents to the tube and repeat magnetic pelleting step.
5. Pipette **50 μ L of Ligation Cleanup Buffer** into each tube to resuspend beads.

Ligation Cleanup Procedure:

! Caution: When vortexing PCR tubes in the subsequent steps, maintain firm pressure on all lids as the contained detergent may allow lids to open.

1. Pipette the entire volume of Ligation Step 2 reaction into the tubes with Ligation Cleanup Beads and Buffer.
2. Mix samples by vortexing.
3. Incubate reactions at room temperature for **5 minutes**.
4. Mix samples by vortexing.
5. Incubate reactions at room temperature for **5 minutes**.
6. Briefly spin down tubes.

7. Place tubes on the magnet for **1 minute or until beads are fully pelleted** against the tube wall.
8. Carefully pipette off and discard supernatant (**100 µL**) without disturbing the beads.
9. Wash beads **two times** with **Ligation Cleanup Buffer**.
 - a. Resuspend beads in **180 µL** Ligation Cleanup Buffer by vortexing, briefly spinning down, and placing back on magnet for **1 minute**.
 - b. Once slurry has cleared, discard supernatant.
10. Wash beads **once** with **ultrapure water**:
 - a. Resuspend beads in **180 µL** of ultrapure water by vortexing, briefly spinning down and placing back on magnet.
 - b. Once slurry has cleared, discard supernatant.
 - c. Take care to ensure that all supernatant has been removed from tubes.
11. Elute DNA from Ligation Cleanup Beads:
 - a. Resuspend ligation cleanup beads in the amount of 5 mM NaOH specified in the panel specific **Product Insert (A)**.
***Note that for automated workflows it is acceptable to add an additional 1 µL 5 mM NaOH to the ligation bead resuspension volume shown above. This extra 1 µL should be left behind when removing the purified DNA in upcoming steps.**
 - b. Transfer reactions to a thermal cycler and initiate an incubation using the following program and guidelines:
 - c. Use a heated lid ($\geq 100^{\circ}\text{C}$).

Ligation Elution incubation conditions

Step	Temperature ($^{\circ}\text{C}$)	Time (minutes)
1	75	10
2	4	Hold

- d. After sample has reached 4°C , briefly spin down and transfer to the magnet.

Step 6: First PCR

1. Thaw the **First PCR Buffer 2.0** and **GSP1** on ice. Briefly vortex and spin down Buffer and primers. Either vortex or pipette mix enzymes then spin down.
2. On ice, make a **First PCR Master Mix** working solution based on the desired number of reactions (including but not exceeding 20% overage) using the table below.
***Note that the PCR Enzyme will also be utilized in Second PCR so do not discard unused enzyme after making First PCR Master Mix.**

First PCR Master Mix	1X (μL)	___X (μL)
First PCR Buffer 2.0 (10029094 or 10029076)	11.2	
PCR Enzyme (10029093 or 10029075)	0.8	
Total volume	12	

- a. Mix by briefly vortexing and spin down.
- b. Keep master mix on ice.
3. To new 0.2mL 8-strip tubes or plate wells add:
 - a. **12 µL** First PCR Master Mix
 - b. GSP1 - Volume specified in the panel specific **Product Insert (B)**
***Note that for automated workflows, the GSPs for the panel(s) in use should be added to their own master mix instead of directly to the sample.**
 - c. Entire volume of purified DNA from Step 5: Ligation 2 Cleanup.
4. Mix by briefly vortexing and spin down.
5. Return tubes to ice.
6. Transfer reactions to a preheated thermal cycler and immediately initiate the program specified in the panel specific **Product Insert (First PCR Reaction)**.
 - a. Use a heated lid ($\geq 100^{\circ}\text{C}$).
 - b. After the program has reached 4°C , briefly spin down reactions and place on ice. It is also acceptable to leave tubes in the thermal cycler at 4°C overnight.

Reaction Cleanup after First PCR

Refer to **Warnings and Precautions** section of this Instructions for Use for guidance on working with Reaction Cleanup beads.

1. Completely resuspend Reaction Cleanup beads by vortexing.
2. Add **0.8X** volume (**32 µL**) of Reaction Cleanup beads to each reaction.
3. Vortex well or pipette 10 times to mix and visually inspect the color of the sample to ensure a homogenous mixture.
4. Incubate for **5 minutes** at room temperature (20°C to 25°C).

5. Briefly spin down tubes.
6. Place tubes on the magnet for **4 minutes or until beads are fully pelleted** against the tube wall.
7. Without disturbing the bead pellet, use a pipette to remove and discard the supernatant. If the pellet becomes dislodged from the magnet and a portion is drawn into the pipette tip, return contents to tube and repeat magnet incubation step.
8. Wash beads **two times** with 70% ethanol while still on the magnet. For each wash:
 - a. Add **180 µL** 70% ethanol
 - b. Incubate for **30 seconds** at room temperature (20°C to 25°C)
 - c. Carefully remove ethanol and discard
9. After the final wash, use a pipette (≤20 µL capacity) to completely remove visible supernatant residue and allow tubes to dry for **3-5 minutes** at room temperature with open lids. **Take care not to over-dry beads** as this will significantly decrease overall recovery (yield) of nucleic acid.
10. Elute DNA by resuspending beads in the volume of 10 mM Tris-HCl, pH 8.0 specified in the panel specific **Product Insert (C)**.
11. Place tubes back on the magnet for **2 minutes**.
12. Transfer the volume of purified eluate in the panel specific **Product Insert (D)** to a new 0.2mL PCR tube and proceed directly to Step 7: Second PCR.

Step 7: Second PCR

1. Thaw the **Second PCR Buffer 2.0**, **GSP2**, and **Liquid P7 Index Plate(s)** (10029147 or 10029150) on ice. Vortex buffer and GSP2 then spin down. Spin down P7 plate before unsealing. It is acceptable to pierce the foil seal with a pipette tip to access the liquid adapters.
2. On ice, make a **Second PCR Master Mix** working solution based on the desired number of reactions (including but not exceeding 20% overage) using the table below.

Second PCR Master Mix	1X (µL)	__X (µL)
Second PCR Buffer 2.0 (10029095 or 10029077)	11.2	
PCR Enzyme (10029093 or 10029075)	0.8	
<i>Total volume</i>	<i>12</i>	

- a. Mix by briefly vortexing and spin down.
 - b. Keep master mix on ice.
3. To new 0.2mL 8-strip tubes or plate wells add:
 - a. **12 µL** Second PCR Master Mix
 - b. **4 µL** respective Liquid P7 Index according to sample index planning. Re-seal P7 plate using a new seal after adding P7 index to respective wells.
 - c. GSP2 - Volume specified in panel specific **Product Insert (E)**.
***Note that for automated workflows, the GSPs for the panel(s) in use should be added to their own master mix instead of directly to the sample.**
 - d. Transfer volume of purified DNA from Step 6: First PCR specified in the panel specific **Product Insert (D)**.
***Important* As this step incorporates the P7 index tag for sample-level tracking, be sure to record which Liquid P7 Index is being used for each sample.**
4. Mix by briefly vortexing and spin down.
5. Return tubes to ice.
6. Transfer reactions to a preheated thermal cycler and immediately initiate the program specified in the panel specific **Product Insert (Second PCR Reaction)**.
 - a. Use a heated lid (≥100°C).
 - b. After the program has reached 4°C, briefly spin down reactions and place on ice. It is also acceptable to leave tubes in the thermal cycler at 4°C overnight.

Reaction Cleanup after Second PCR

Refer to **Warnings and Precautions** section of this Instructions for Use for guidance on working with Reaction Cleanup beads.

1. Completely resuspend Reaction Cleanup beads by vortexing.
2. Add **0.8X** volume (**32 µL**) of Reaction Cleanup beads to each Second PCR reaction.
3. Vortex well or pipette 10 times to mix and visually inspect the color of the sample to ensure a homogenous mixture.
4. Incubate for **5 minutes** at room temperature (20°C to 25°C).
5. Briefly spin down tubes.
6. Place tubes on the magnet for **4 minutes or until beads are fully pelleted** against the tube wall.

7. Without disturbing the bead pellet, use a pipette to remove and discard the supernatant. If the pellet becomes dislodged from the magnet and a portion is drawn into the pipette tip, return contents to tube and repeat magnet incubation step.
 8. Wash beads **two times** with 70% ethanol while still on the magnet. For each wash:
 - a. Add **180 µL** 70% ethanol.
 - b. Incubate for **30 seconds** at room temperature (20°C to 25°C).
 - c. Carefully remove ethanol and discard.
 9. After the final wash, use a pipette (≤20 µL capacity) to completely remove visible supernatant residue and allow tubes to dry for **3-5 minutes** at room temperature with open lids. **Take care not to over-dry beads** as this will significantly decrease overall recovery (yield) of nucleic acid.
 10. Elute DNA by resuspending beads in **20 µL** 10 mM Tris-HCl, pH 8.0.
 11. Place tubes back on the magnet for **2 minutes**.
 12. Transfer **18 µL** of the purified solution to a new 0.2mL PCR tube. Be sure to avoid transferring beads to the fresh tube.
 13. Stop or proceed directly to Step 8: Quantify Libraries.
-  **Safe stopping point:** It is okay to stop and store the reactions at -30°C to -10°C.

Step 8: Quantify Libraries

1. Quantify the concentration of each library by qPCR following the instructions provided with the KAPA Universal Library Quantification Kit.
-  **Note:** Archer libraries are very concentrated and may need to be serially diluted to 1:100,000 for quantification with KAPA via qPCR.
2. The Loading Calculator spreadsheet can be used for quantifying, normalizing, and sequencing calculations.
 3. The recommended size-adjustment correction factor for calculating the concentration of a library is dependent on the input type and quality. The estimated fragment size values below account for adapter size (157 bp) plus the average insert fragment length. Use an estimated fragment size of 357 bp.
 4. Use each library's quantity to guide dilution of an aliquot to enable combining all normalized libraries to create an equimolar pool (e.g, 10nM).
-  **Note:** Quantification with a fluorometric or electrophoresis-based method is not recommended.

Quality Control

This product was manufactured and released in accordance with the Integrated DNA Technologies, BV, ISO 13485 Quality Management System. The product has met Integrated DNA Technologies quality requirements.

Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and eye protection. For more information, please consult the specific safety data sheet (SDS) by searching on our website at <https://integrateddna-sds.thewerco.com/>.

Disposal

Dispose of reagent and packaging in accordance with local regulations.

Glossary of Symbols

A glossary of symbol definitions is available for reference at www.idtdna.com.

Symbol	Description
	English language
	In vitro diagnostic medical device
	Conformité Européene (CE) Mark
	Authorized representative in the European Community/European Union
	Manufacturer's catalogue number.
	Manufacturer's batch code.
	The date after which the product is not to be used.
	Indicates the manufacturer.
	The date when the product was manufactured.

	The total number of tests that can be performed with the product.
	The temperature limits to which the product can be safely exposed.
	User to consult the instructions for use.
	Reminder symbols call attention to minor details that may be easily overlooked and compromise the procedure resulting in decreased assay performance.
	Caution symbols denote critical steps in the procedure where risk of protocol failure or damage to the product itself could occur if not carefully observed.
	Stop symbols indicate where this procedure may be safely suspended and resumed at a later time without risk of compromised assay performance. Make note of these steps and plan your workflow accordingly.
	French Language
	Italian Language
	Spanish Language

Version History

Version Number	Version Date	Description of Change
01	04/2025	New document
01	11/2025	Converted to new document type (PRD-DOC-391 to RA-DOC-625). Condensed content, table formats, and margins for print version. Added to definitions table. Added foreign language translations (French, Spanish, and Italian).

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