


xGen™ cfDNA & FFPE DNA Library preparation for sequencing using Ultima Genomics UG 100

The method presented here is provided by IDT. This method uses the xGen cfDNA & FFPE DNA Library Prep v2 MC Kit and xGen Indexing Primers for Ultima with Ultima Genomics' UG Library Amplification mix. This method can be used as a starting point for creating next generation sequencing (NGS) libraries intended for downstream UG 100 sequencing, and can be used in similar experiments. However, it may not be fully optimized for your application. IDT does not guarantee these methods, and application scientists at IDT can only provide general guidance with limited troubleshooting support. For more information and support for Ultima products and sequencing on the UG 100 System, see Ultima Genomics' website www.ultimagenomics.com or email support@ultimagenomics.com

Overview

The xGen cfDNA & FFPE DNA Library Prep v2 MC Kit is designed specifically for 1–250 ng of degraded samples, such as cell-free DNA (cfDNA) or DNA extracted from formalin-fixed, paraffin-embedded (FFPE) samples. The technology utilized for this library construction workflow enables high conversion of degraded or damaged samples into sequencer-ready libraries without the need for adapter titration. The kit also allows for independent tagging of the top and bottom strands of the library, using our proprietary ligation strategy with fixed UMI sequences, which enable bioinformatic error correction. Combining higher conversion rate with error correction better enables identification of ultra-low frequency variants.

This demonstrated protocol describes the workflow for preparing libraries for whole genome sequencing (WGS) on the UG 100 sequencing platform from Ultima Genomics. For sequencing, the UG 100 system requires libraries containing a UG 100 specific Index Primer and an UG Universal Primer sequence containing an overhang.

 **Important:** Therefore, the UG Library Amplification Kit and xGen Indexing Primers for Ultima must be used with the xGen cfDNA & FFPE DNA Library Prep v2 MC Kit to create compatible libraries for the UG 100 sequencing platform.

There are five major steps to perform in this protocol:

- **End repair.** Converts cfDNA, or sheared input DNA, into blunt-ended DNA ready for ligation.
- **Ligation 1.** Ligation 1 Enzyme catalyzes the single-stranded addition of the Ligation 1 Adapter to the 3' ends of the insert DNA.
- **Ligation 2.** Ligation 2 Adapter acts as a primer to gap-fill the bases complementary to the UMI, followed by ligation to the 5' end of the DNA insert to create a fully double-stranded product.
- **PCR amplification.** UG Library Amplification Kit and xGen Indexing Primers for Ultima are used to perform indexing PCR for sequencing.
- **Overhang generation.** Library end is nicked to create an overhang to allow for binding to sequencing beads.

> SEE WHAT MORE WE CAN DO FOR YOU AT WWW.IDTDNA.COM.

| | | | |
|---|--|-------------|--|
| 1 | End repair | | Total time: 35 min |
| 2 | End repair cleanup | 2.5X AMPure | Total time: 30 min |
| 3 | Ligation 1 | | Total time: 40 min |
| 4 | Ligation 2 | | Total time: 45 min |
| 5 | Ligation 2 cleanup | 2.5X | Total time: 45 min |
| 6 | PCR amplification using UG specific reagents | | Total time: 15–30 min (variable depending on input) |
| 7 | PCR | 1.2X AMPure | Total time: 40 min |
| 8 | Overhang generation | | Total time: 40 min |
| 9 | Overhang generation cleanup | 1.2X AMPure | Total time: 40 min |

Consumables and equipment

Consumables from IDT

| Item | Catalog # |
|---|-------------|
| xGen cfDNA & FFPE DNA Library Prep v2 MC Kit | 10010207 |
| Indexing Primers for Ultima (choose one): | |
| • xGen Indexing Primers for Ultima, 16 rxn | 10016991 |
| • xGen Indexing Primers for Ultima P1, 96 rxn | 10016992 |
| • xGen Indexing Primers for Ultima P2, 96 rxn | 10016993 |
| • xGen Indexing Primers for Ultima P3, 96 rxn | 10016994 |
| • xGen Indexing Primers for Ultima P4, 96 rxn | 10016995 |
| IDTE pH 8.0 (1X TE Solution), 300 mL | 11-05-01-13 |

Consumables—Other suppliers

| Item | Supplier | Catalog # |
|--|--------------------------|------------------|
| UG Library Amplification Kit v3.0 | Ultima Genomics | 20020007 |
| Buffer EB (10 mM Tris-HCl, pH 8.5), 250 mL | QIAGEN | 19086 |
| Absolute ethanol (200 proof) | Various | Varies |
| Purification beads | | |
| Agencourt® AMPure® XP-PCR Purification Beads 5 mL (22 rxn) or 60 mL (270 rxn) | Beckman Coulter | A63880 or A63881 |
| Digital electrophoresis chips and associated reagents (choose one): | | |
| • Experion® DNA 1K Analysis Kit | Bio-Rad | 700-7107 |
| • High Sensitivity DNA Kit | Agilent | 5067-4626 |
| • High Sensitivity D1000 ScreenTape® | Agilent | 5067-5584 |
| twin.tec® 96 Well LoBind® PCR Plates | Eppendorf | 0030129504 |
| DNA LoBind® Tubes, 1.5 mL | Eppendorf | 022431021 |
| Filtered low bind pipette tips | Various | Varies |
| Qubit™ dsDNA HS Assay Kit | Thermo Fisher Scientific | 032851 or 032854 |

Equipment

| Item | Supplier | Catalog # |
|--|--------------------------|-----------|
| Microcentrifuge | Various | Varies |
| Thermal cycler | Various | Varies |
| Qubit 4 Fluorometer, or similar DNA quality analyzer | Thermo Fisher Scientific | 033226 |
| Magnet options (choose one): | | |
| • Permagen 96-well side pull magnet | TMO | NC1568572 |
| • Magnum® EX Universal Magnet Plate | Alpaqua | A000380 |
| • Magnetic Stand-96, or similar magnetic stand | Thermo Fisher Scientific | AM10027 |

Guidelines

Reagent storage and handling


Always store the xGen cfDNA & FFPE DNA Library Prep v2 MC Kit reagents at -20°C , except for PEG/NaCl, which can be stored either at room temperature, or at -20°C .


Notes:

- The enzymes provided in this kit are temperature-sensitive and appropriate care should be taken during storage and handling.
- The xGen 2x HiFi PCR Mix may be in a liquid state after storage in -20°C condition; this is expected and will not impact product integrity.

For all non-enzyme reagents, thaw on ice, then briefly vortex to mix well. Remove enzyme tubes from -20°C storage and place on ice just before use. Briefly centrifuge the tubes in a microcentrifuge to collect contents before opening.


To create master mixes, scale reagent volumes as appropriate, using 10% excess volume to compensate for pipetting loss. Add reagents in the order listed when preparing master mixes, then pulse-vortex to mix and briefly centrifuge.

 **Note:** You may observe precipitation in the End Repair Buffer. Continue vortexing until no precipitate can be observed.

 **Important:** Use extra caution when handling Ligation 1 and Ligation 2 adapter tubes. Never handle Ligation 2 Adapter prior to, or during Ligation 1 Master Mix setup and handling. Trace contamination of Ligation 2 Adapter into Ligation 1 Adapter has been shown to induce adapter-dimer formation.

DNA input considerations

This kit works with a wide range of DNA inputs, ranging from 1–250 ng. Input DNA should be stored in IDTE buffer (10 mM Tris-Cl, pH 8.0; 0.1 mM EDTA) or Buffer EB (QIAGEN, 10 mM Tris HCl, pH 8.5).

 **Important:** Input quantities recommended in this protocol refer to the total DNA quantified after fragmentation.

DNA should be an appropriate size before library construction. For genomic DNA, or DNA derived from FFPE samples, using a Covaris® or similar DNA shearing instrument can create fragments with an average insert size of 150–300 base pairs. As cfDNA typically has an average size of 160 base pairs, no further fragmentation is required.

For FFPE samples, use standard quality control methods, such as Q-ratio with qPCR or the DNA Integrity Number (DIN) using size distribution (i.e., Bioanalyzer instrument (Agilent) or similar DNA quality analyzer).

These methods can help you choose the appropriate number of PCR cycles for your DNA sample. For cfDNA, we suggest assessing the size distribution with electrophoresis. Note, if large molecular weight DNA is present, an additional cleanup may be necessary to remove genomic DNA contamination. However, this additional cleanup can reduce sample complexity and mass.

 **Note:** For more information on avoiding cross contamination, and size selection during clean-ups see the full [xGen cfDNA & FFPE DNA Library Prep v2 MC protocol](#).

Protocol

End repair

Before starting the protocol, ensure that the AMPure and PEG/NaCl reagents are at room temperature (20–25°C). Also, prepare a fresh 80% ethanol solution in nuclease-free water.

1. Add 50 µL of each sample into a low-bind PCR plate that resists nucleic acid adsorption.



Note: If sample volume <50 µL, use IDTE buffer (10 mM Tris-Cl, pH 8.0; 0.1 mM EDTA) or Buffer EB (10 mM Tris-HCl, pH 8.5) to bring the volume up to 50 µL.

2. For each sample, make the following End Repair Master Mix:

| End Repair Master Mix | |
|-----------------------|--------------------------|
| Component | Volume per reaction (µL) |
| End Repair Buffer | 6 |
| End Repair Enzyme | 3 |
| Total volume: | 9 |



Tip: If there is precipitate in the End Repair Buffer, vortex until the precipitate becomes clear in solution. The resulting Master Mix is viscous and requires careful pipetting.

3. Pulse-vortex the master mix for 10 sec, then briefly centrifuge. Keep the master mix on ice.
4. Add 9 µL of End Repair Master Mix to each well. Using a pipette set to 40 µL, pipette 10 times to mix, then seal the plate. Alternatively, seal the plate and vortex mix for minimum of 10 seconds.
5. Run the following thermal cycler program:

| End Repair program | | |
|--------------------|-------------------|--------|
| Step | Temperature* (°C) | Time |
| End repair | 20 | 30 min |
| Hold | 4 | Hold |

* Set the lid temperature to OFF, or to 40°C.

While the end repair program runs, make the Ligation 1 Master Mix in preparation for the End repair cleanup steps.

| Ligation 1 Master Mix | |
|-----------------------|--------------------------|
| Component | Volume per reaction (µL) |
| Ligation 1 Buffer | 25 |
| Ligation 1 Adapter | 2 |
| Ligation 1 Enzyme | 3 |
| Total volume: | 30 |

6. Pulse vortex the master mix for 10 seconds, then briefly centrifuge. Keep the Master Mix on ice until ready to use.




Important: Use extra caution when handling Ligation 1 and Ligation 2 Adapter tubes: never handle Ligation 2 Adapter before or during Ligation 1 Master Mix setup and handling. Trace contamination of Ligation 2 Adapter into Ligation 1 Adapter has been shown to induce adapter-dimer formation.


7. After the End Repair program reaches 4°C, proceed immediately to [End repair cleanup](#).


End repair cleanup

 **Note:** Before starting cleanup, make sure the Ligation 1 Master Mix has already been prepared.

1. Thoroughly resuspend AMPure XP beads before use, then add 147.5 μ L of AMPure beads (2.5X volume) to each well and pipette 10 times to thoroughly mix.
2. Incubate the plate at room temperature for 10 minutes.
3. Place the plate on a magnet and wait for the liquid to clear completely or for at least 2 minutes.


 **Important:** If the solution is not clear after 2 minutes, keep the plate on the magnet until solution clears. If beads are discarded, sequencing library yield and quality may suffer.
4. Remove and discard the cleared supernatant making sure not to remove beads. With a clean P20 pipette tip, remove and discard any trace amount of supernatant that remains.

 **Important:** Be careful not to remove any beads.
5. Keeping the plate on the magnet, add 160 μ L of 80% ethanol and incubate for 30 seconds.
6. Remove and discard the supernatant.
7. Use a P20 pipette tip to remove any residual ethanol.
8. Dry the beads at room temperature for 1–3 minutes.

 **Important:** Never allow the beads to dry for >5 minutes; over-dry beads will decrease final library yield.
9. Proceed immediately to [Ligation 1](#).


Ligation 1

1. Remove the plate from the magnet, then add 30 μ L Ligation 1 Master Mix.
2. Pipette mix a minimum of 10 times, then seal the plate.

 **Important:** Make sure the samples are thoroughly mixed and that the beads are fully resuspended before proceeding.
3. Run the following thermal cycler program:

| Ligation 1 program | | |
|--------------------|-------------------|--------|
| Step | Temperature* (°C) | Time |
| Ligation | 20 | 15 min |
| Inactivation | 65 | 15 min |
| Hold | 4 | Hold |

*If possible, set lid temperature to 70°C. If lid cannot be programmed, set the thermal cycler heat block to 105°C.



 **Safe Stop:** The plate can temporarily remain at 4°C (no more than 2 hours). It is normal for beads to settle during this reaction.

4. Proceed to [Ligation 2](#).

Ligation 2

1. For each sample, prepare the Ligation 2 Master Mix.

| Ligation 2 Master Mix | |
|-----------------------|--------------------------|
| Component | Volume per reaction (μL) |
| Ligation 2 Buffer | 4.5 |
| Ligation 2 Adapter | 4 |
| Ligation 2 Enzyme A | 0.5 |
| Ligation 2 Enzyme B | 1 |
| Total volume: | 10 |





2. Pulse-vortex the master mix for 10 seconds, then briefly centrifuge. Keep the master mix on ice until ready to use.
3. Add 10 μL of the Ligation 2 Master Mix to each well.
4. Using a pipette set to 35 μL, pipette 10 times to mix, then seal the plate.
 -  **Important:** Ensure the samples are thoroughly mixed, and make sure the beads are fully resuspended before proceeding.
 -  **Note:** If necessary, briefly centrifuge to collect contents to the bottom of the wells.
5. Run the following thermal cycler program:

| Ligation 2 program | | |
|--------------------|-------------------|--------|
| Step | Temperature* (°C) | Time |
| Ligation | 65 | 30 min |
| Hold | 4 | Hold |

*If possible, set lid temperature to 70°C. If lid cannot be programmed, set the thermal cycler heat block to 105°C.

6. After thermal cycler program is completed, proceed immediately to [Ligation 2 cleanup](#).

Ligation 2 cleanup

1. Add 100 μL of PEG/NaCl (2.5X volume) to each well, then pipette 10 times to mix.
2. Incubate the plate at room temperature for 10 minutes.
3. Place the plate on a magnet and wait for the liquid to clear completely or at least for 2 minutes.
 -  **Important:** If solution is not clear after 2 minutes, keep the plate on the magnet until solution clears. If beads are discarded, sequencing library yield and quality may suffer.
4. Remove and discard the cleared supernatant making sure not to remove beads. Use a clean P20 pipette tip to remove and discard any trace amount of supernatant that remains.
 -  **Important:** Be careful not to remove any beads.
5. Keeping the plate on the magnet, add 160 μL of 80% ethanol and incubate for 30 seconds.
6. Remove and discard the supernatant.
7. Use a P20 pipette tip to remove any residual ethanol.
 -  **Important:** Make sure all the ethanol has been removed before proceeding.
8. Dry the beads at room temperature for 1–3 minutes.
 -  **Important:** Never allow the beads to dry for >5 minutes; over-dry beads will decrease final library yield.

9. Remove the plate from the magnet, then add 20 μ L of Buffer EB.
10. Seal the plate, then gently vortex (use 70% vortex capacity) to resuspend beads.
11. Allow the plate to incubate at room temperature for 5 minutes to elute DNA off the beads.
12. Place the plate on a magnet and wait for the beads to be cleared from the liquid (approximately 1–2 minutes).



Note: Depending on the strength of your magnet, you may need to wait longer.

13. Carefully transfer 20 μ L of the cleared liquid containing the eluted DNA into a new well. Proceed to **PCR amplification** or pause here.



Safe Stop: The plate can be stored at -20°C overnight.

PCR amplification using UG specific reagents

Important: UG Library Amplification Kit and xGen Indexing Primers for Ultima must be used during PCR amplification. Do not use the amplification mix provided in the xGen cfDNA & FFPE DNA Library Prep v2 MC Kit.

Note: Sample index barcode is introduced during PCR; double check that a unique primer pair is used for each sample.

1. Add 5 μ L xGen Indexing Primers for Ultima to each well.
2. Add 25 μ L of UG Master Mix to each well, then pipette 10 times to mix.
3. Seal the plate, then briefly centrifuge.
4. Run the following thermal cycler program:



Note: Only use standard cycling conditions for the PCR program. Fast cycling has been shown to negatively impact library yields.

| PCR program | | | |
|--------------------------|-------------------------------------|--------|---|
| Step | Temperature* ($^{\circ}\text{C}$) | Time | Cycles |
| Polymerase activation | 98 | 45 sec | Varies based on sample input (see Table 1) |
| Denature, anneal, extend | 98 | 15 sec | |
| | 60 | 30 sec | |
| | 72 | 30 sec | |
| Final extension | 72 | 1 min | |
| Hold | 4 | Hold | |

*Set lid temperature to 105°C . If the lid cannot be programmed, set the thermal cycler block to 105°C .

Table 1. Recommended PCR cycling parameters to yield >500 ng.


| Input mass (ng) | Number of cycles | |
|-----------------|------------------|-------|
| | gDNA or dsDNA | FFPE* |
| 1 | 11–13 | 14–16 |
| 10 | 9–11 | 11–13 |
| 25 | 7–9 | 9–11 |
| 100 | 5–7 | 7–9 |
| 250 | 4–6 | 6–8 |


* For lower-quality FFPE samples, we recommend using the higher number of recommended PCR cycles from [Table 1](#). For very low quality FFPE (DIN 1–2), you may not obtain 500 ng, but we do not recommend increasing the number of PCR cycles.


5. After the program completes, proceed to **PCR cleanup**.


PCR cleanup


1. Add 60 μ L of AMPure beads (1.2X volume), or similar bead-based DNA purification product, to each well, then pipette 10 times to thoroughly mix.
2. Incubate the plate at room temperature for 5 minutes.
3. Place the plate on a magnet and wait for the liquid to clear completely or at least for 2 minutes.

 **Important:** If solution is not clear after 2 minutes, keep the plate on the magnet until solution clears. If beads are discarded, sequencing library yield and quality may suffer.
4. Remove and discard the cleared supernatant; make sure not to remove any beads.

 **Important:** Be careful not to remove any beads.
5. Keeping the plate on the magnet, add 200 μ L of 80% ethanol, then incubate for 30 seconds.
6. Remove and discard the supernatant.
7. Repeat the wash step for a total of two washes.
8. Use a P20 pipette tip to remove any residual ethanol.

 **Important:** Make sure all the ethanol has been removed before proceeding.
9. Dry the beads at room temperature for 1–3 minutes.

 **Important:** Never allow the beads to dry for >5 minutes; over-dry beads will decrease final library yield.
10. Remove the plate from magnet, then add 31 μ L of Buffer EB.
11. Seal the plate and gently vortex (use 70% vortex capacity) to resuspend beads.
12. Allow the plate to incubate at room temperature for 5 minutes to elute DNA off beads. Then place the plate on a magnet and wait for the liquid to clear completely for 1–2 minutes.

 **Note:** Depending on the strength of your magnet, you may need to wait longer.
13. Carefully transfer 30 μ L of eluted DNA into a new well.

Overhang generation using UG specific reagents

1. For each sample, prepare the Nicking Master Mix.

| Nicking Master Mix | |
|----------------------|--------------------------------|
| Component | Volume per reaction (μ L) |
| • Nicking Mix Buffer | 4 |
| • UG Nicking Mix | 6 |
| Total volume: | 10 |







2. Pulse vortex the master mix for 10 seconds, then briefly centrifuge. Keep the master mix on ice until ready to use.
3. Add 10 μ L of the Nicking Master Mix to 30 μ L of purified amplified library.
4. Seal the plate, mix thoroughly by vortex, and briefly centrifuge.
5. Run the following thermal cycler program:

| Step | Temperature* ($^{\circ}$ C) | Time |
|------|------------------------------|------------|
| 1 | 37 | 30 minutes |
| 2 | 4 | Hold |

* Lid temperature needs to be set to 40 $^{\circ}$ C or off. If lid cannot be programmed, set to 105 $^{\circ}$ C

6. After thermal cycler program is completed, proceed immediately to [Overhang generation clean-up](#).

Overhang generation cleanup

1. Add 48 µL of AMPure beads (1.2X volume), or similar bead-based DNA purification product, to each well, then vortex to mix, pulse spin to bring all liquid to the bottom of well (do not pellet beads).
2. Incubate the plate at room temperature for 5 minutes.
3. Place the plate on a magnet and wait for the liquid to clear completely or at least for 2 minutes.
 -  **Important:** If solution is not clear after 2 minutes, keep the plate on the magnet until solution clears. If beads are discarded, sequencing library yield and quality may suffer.
4. Remove and discard the cleared supernatant; make sure not to remove any beads.
 -  **Important:** Be careful not to remove any beads.
5. Keeping the plate on the magnet, add 200 µL of 80% ethanol, then incubate for 30 seconds.
6. Remove and discard the supernatant.
7. Repeat wash steps for a total of 2 washes.
8. Use a P20 pipette tip to remove any residual ethanol.
 -  **Important:** Make sure all the ethanol has been removed before proceeding.
9. Dry the beads at room temperature for 1–4 minutes.
 -  **Important:** Never allow the beads to dry for >5 minutes; over-dry beads will decrease final library yield.
10. Remove the plate from magnet, then add 31 µL of Buffer EB.
11. Seal the plate and gently vortex (use 70% vortex capacity) to resuspend beads.
12. Allow the plate to incubate at room temperature for 5 minutes to elute DNA off beads. Then place the plate on a magnet and wait for the liquid to clear completely for 1–2 minutes.
 -  **Note:** Depending on the strength of your magnet, you may need to wait longer.
13. Carefully transfer 30 µL of eluted DNA into a new well.
 -  **Safe Stop:** Libraries can be stored at –20°C overnight.

The library is now ready for quantification and quality check. For information on library QC, pooling and recommended loading concentrations please consult the Ultima Genomics Library QC and Pooling Reference Guide (P00045).

Library performance

DNA was extracted from FFPE samples (SeraSeq) using the QIAGEN QiAmp DNA extraction kit and DIN scores were determined using the TapeStation. Example libraries were prepared from 25 ng input of DNA using the xGen cfDNA & FFPE DNA Library Prep v2 MC Kit with the UG Library Amplification Kit and custom IDT indexing primers equivalent to the xGen Indexing Primers for Ultima P1. In [Figure 1](#), a TapeStation trace (Agilent) shows a PCR-amplified library generated from a low-quality FFPE sample of DIN 2. Even with the low DIN score a quality library was generated without any evidence of adapter dimer. [Table 2](#) shows the quality sequencing metrics obtained from sequencing 12 xGen cfDNA & FFPE libraries on the UG 100 sequence. The libraries had a high percentage of reads aligned to the reference genome, and consistent results across samples with varying DIN levels. The data shows low indel, mismatch, and chimera rates and high Q20 indicating that quality sequencing reads were generated on the UG 100 system. Table represents the mean of the three replicates.




Figure 1. Representative TapeStation trace for xGen cfDNA & FFPE DNA library compatible with UG 100 sequencer. The FFPE sample had a DIN score of 2. The resulting library showed a peak height of ~325 bp; this is the expected library size after adding the adapters and index sequences.

Table 2. Quality metrics from IDT libraries sequenced on the UG 100 system.

| | Aligned reads | Indel error rate | Mismatch error rate | Chimera | Bases>Q20 |
|------------------|---------------|------------------|---------------------|---------|-----------|
| Sample 1 DIN 2 | 98.40% | 0.30% | 0.40% | 1.10% | 96.40% |
| Sample 2 DIN 5.8 | 98.90% | 0.30% | 0.30% | 0.50% | 95.80% |
| Sample 3 DIN 5.8 | 98.90% | 0.30% | 0.30% | 0.20% | 95.30% |

Fixed UMI sequences

The Ligation 1 Adapter contains 32 optimized, fixed UMI sequences that are 8 base pairs in length.

| Adapter name | UMI sequence | Adapter name | UMI sequence |
|--------------|--------------|--------------|--------------|
| UMI_1 | GAGACGAT | UMI_17 | GCACAAC |
| UMI_2 | TTCCAAGG | UMI_18 | GCGTCATT |
| UMI_3 | CGCATGAT | UMI_19 | GAAGGAAG |
| UMI_4 | ACGGAACA | UMI_20 | ACTGAGGT |
| UMI_5 | CGGCTAAT | UMI_21 | TGAAGACG |
| UMI_6 | GCTATCCT | UMI_22 | GTTACGCA |
| UMI_7 | TGGACTCT | UMI_23 | AGCGTGTT |
| UMI_8 | ATCCAGAG | UMI_24 | GATCGAGT |
| UMI_9 | CTTAGGAC | UMI_25 | TTGCGAAG |
| UMI_10 | GTGCCATA | UMI_26 | CTGTTGAC |
| UMI_11 | TCGCTGTT | UMI_27 | GATGTGTG |
| UMI_12 | TTCGTTGG | UMI_28 | ACGTTTCA |
| UMI_13 | AAGCACTG | UMI_29 | TTGCAGAC |
| UMI_14 | GTCGAAGA | UMI_30 | CAATGTGG |
| UMI_15 | ACCACGAT | UMI_31 | ACGACTTG |
| UMI_16 | GATTACCG | UMI_32 | ACTAGGAG |

xGen™ cfDNA & FFPE DNA Library preparation for sequencing using Ultima Genomics UG 100

For more information, go to: www.idtdna.com/ContactUs

For more than 30 years, IDT's innovative tools and solutions for genomics applications have been driving advances that inspire scientists to dream big and achieve their next breakthroughs. IDT develops, manufactures, and markets nucleic acid products that support the life sciences industry in the areas of academic and commercial research, agriculture, medical diagnostics, and pharmaceutical development. We have a global reach with personalized customer service.

> SEE WHAT MORE WE CAN DO FOR YOU AT WWW.IDTDNA.COM.

For Research Use Only. Not for use in diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2024 Integrated DNA Technologies, Inc. All rights reserved. Trademarks contained herein are the property of Integrated DNA Technologies, Inc. or their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks.
Doc ID: RUO23-2555_001 01/24