# **B** Element Biosciences



# Quick Guide for the Trinity Workflow xGen<sup>™</sup> Exome Hybridization & Trinity Run Setup

# Introduction

Intended for the experienced user, this quick guide describes the xGen<sup>™</sup> Exome hybridization protocol and Trinity sequencing run setup on an AVITI<sup>™</sup> System.

The protocol supports up to a 24-plex pooled hybridization reaction in one capture reaction.

# Prepare Libraries for Hybridization

- 1. Gather the following components for hybridization:
  - » Elevate linear library, prepared using the xGen<sup>™</sup> DNA Library Prep EZ UNI Kit
  - » xGen<sup>™</sup> Hyb Capture Kit Trinity<sup>™</sup> for Element
    - ∘ xGen<sup>™</sup> Human Cot DNA
    - xGen<sup>™</sup> 2X Hybridization Buffer
    - xGen<sup>™</sup> Hybridization Buffer Enhancer
    - ∘ xGen<sup>™</sup> Exome Hyb Panel
  - » Trinity Binding Reagent
- 2. Thaw reagents on ice and mix thoroughly.
- 3. Add the total library volume and the following reagents to a 1.5 ml LoBind tube or a 96-well plate. Vortex to mix and briefly centrifuge.

Component	Volume
Up to 24 indexed samples	5–6 μg total
xGen™ Human Cot DNA	5 µl
Trinity Binding Reagent	5 µl

- 4. Dry the indexed pool using a SpeedVac at low or no heat (less than 40°C).
- 5. If not proceeding to hybridization, store the dried-down library pool at -25°C to -15°C up to 1 week.

# Perform Hybridization

 Inspect the xGen<sup>™</sup> 2X Hybridization Buffer for salt crystals. If present, heat the tube at 65°C and shake intermittently until the buffer is solubilized.

#### Prepare the hybridization master mix in a 1.5 ml tube. Multiply volumes by the number of samples and add 10% overage.

Component	Volume
xGen™ 2X Hybridization Buffer	8.5 μl
xGen™ Hybridization Buffer Enhancer	2.7 μl
xGen™ Exome Hyb Panel	4 µl
Water	1.8 µl
Total	17 µl

- 3. Vortex or pipette to mix the hybridization master mix solution and then briefly centrifuge.
- 4. Transfer 17  $\mu$ l hybridization master mix to each tube or well containing dried libraries. Pipette to mix.
- 5. Incubate at room temperature for 5–10 minutes.
- 6. If using tubes up to this step, transfer the full volume to a 96-well PCR plate.
- 7. Cap the tube or seal the plate tightly to avoid evaporation. Vortex and then briefly centrifuge.
- 8. Run the following thermal cycler program to incubate.

Temperature	Time
Lid set to 100°C	
95°C	30 seconds
65°C	16 hours
65°C	Hold

## **Thaw Reagents**

1. Thaw the Trinity sequencing cartridge. Protect from light.

Cartridge	Water Bath	Refrigerator
2 x 75	90 minutes	8 hours
2 x 150	2.5 hours	24 hours

- 2. Make sure reagents are *fully* thawed.
- 3. Set aside at room temperature or keep at 2°C to 8°C.

#### xGen<sup>™</sup> Exome Hybridization & Trinity Run Setup

# Initiate a Sequencing Run

- 1. On the Home screen, select New Run.
- 2. Select Sequencing.
- 3. Select the side for sequencing:Side A, Both, or Side B.
- 4. For chemistry type, select **Trinity**, and then select **Next**.
- For a Manual Run, proceed to <u>Define Run Parameters</u>.
  For a Planned Run, select the run and storage connection, and then select Next. Proceed to <u>Inspect and Mix</u> <u>Reagents</u>.

### **Define Run Parameters**

- 1. In the Run Name field, enter a unique name.
- 2. If applicable, select **Browse** and import the run manifest.
- 3. Complete the Description and Storage fields as applicable.
- 4. Select a Trinity Sequencing Kit.
- 5. Select the panel **xGen Exome Kit for Trinity**.
- 6. Enter the number of cycles, and then select Next.

### Inspect and Mix Reagents

- 1. Gently invert the cartridge *10 times*.
- 2. Tap the base on the benchtop.
- 3. Place into a cartridge basket and lock the clips.

### **Prepare Sequencing Solution**

- 1. Gather the following components:
  - » Trinity Sequencing Reagent
  - » Library Loading Buffer
- 2. Remove the hybridization reaction from the thermal cycler and briefly centrifuge.
- 3. *Immediately* add 183 μl Library Loading Buffer to dilute each xGen<sup>™</sup> hybridization reaction. Pipette gently to mix.
- 4. Prepare the sequencing solution in a 5 ml tube. Pipette gently to mix.

Component	Volume
Library Loading Buffer	2038 µl
Trinity Sequencing Reagent	72 µl
Diluted hybridization reaction	90 µl
Total	2200 μl

## Add Sequencing Solution to Cartridge

- 1. Using a 1 ml pipette tip, pierce the Library well.
- 2. Transfer 2200 µl sequencing solution to the Library well.

## **Confirm Reagent Preparation**

- 1. Select the **Invert cartridge** checkbox.
- 2. Select the **Insert into basket** checkbox.
- 3. Select the Load hybed reaction checkbox. Select Next.

#### Load Reagents and Buffer

- 1. Open the reagent bay door and remove any materials.
- 2. Slide the basket into the reagent bay.
- 3. Slide the buffer bottle into the reagent bay until it stops.
- 4. Close the reagent bay door, and then select Next.

#### **Empty Waste and Prime Reagents**

- 1. Open the waste bay door, remove the waste bottle, and close the transport cap.
- 2. Open the transport and vent caps and empty the waste.
- 3. Close the vent cap and reload the waste bottle.
- 4. Select **Next** to *automatically* start priming.
- 5. Bring a new Trinity flow cell to room temperature in the package.
- 6. When priming is complete, select Next.

#### Load the Flow Cell

- 1. Remove the used flow cell from the nest.
- 2. Unpackage the Trinity flow cell and load it onto the nest.
- 3. Select Close Nest, and then select Next.

### Review and Start the Run

- 1. Review the run, and then select **Run**.
- 2. Monitor run metrics as they appear onscreen.

© 2024 Element Biosciences, Inc. and Integrated DNA Technologies, Inc. Element Biosciences, AVITI, Trinity, and the Element Biosciences logo are trademarks of Element Biosciences, Inc. xGen is a trademark of Integrated DNA Technologies, Inc. and registered in the US. Information in this document is subject to change without notice.

#### xGen<sup>™</sup> Exome Hybridization & Trinity Run Setup