

Product Insert FusionPlexTM Lung Panel

FusionPlex Lung

Description

The FusionPlex Lung is a balanced pool of gene-specific primer (GSP) oligonucleotides that is optimized for use with FusionPlex reagents and molecular barcode (MBC) adapters to produce targeted NGS libraries. This product insert should be used in conjunction with FusionPlex protocol for Illumina® (RA-DOC-047) or FusionPlex protocol for Ion Torrent™ (RA-DOC-048).

FusionPlex Lung contains **163** GSPs targeting **14** genes commonly mutated in lung cancer.

Description	Part number	Storage
FusionPlex Lung GSP1, 8 reactions	SA0166081	
FusionPlex Lung GSP2, 8 reactions	SA0166082	−20°C ± 10°C
10X VCP Primer Mix	SA0126	

Recommended PCR cycling

	Step	Temperature (°C)	Time	Cycles
First PCR reaction	1	95	3 min	1
	2	95	30 sec	
	3	65	5 min (100% ramp rate)	15
	4	72	3 min	1
	5	4	Hold	1

2425 55th Street, Boulder, CO 80301 | www.archerdx.com | archer-tech@idtdna.com



Product Insert FusionPlexTM Lung Panel

Recommended PCR cycling (cont.)

	Step	Temperature (°C)	Time	Cycles	
Second PCR reaction	1	95	3 min	1	
	2	95	30 sec		
	3	65	5 min (100% ramp rate)	20†	
	4	72	3 min	1	
	5	4	Hold	1	

[†]The number of PCR2 cycles may be decreased if you regularly experience library yields greater than 200 nM.

Recommended reads and multiplexing

FusionPlex Lung libraries should be sequenced to a minimum of **0.5M reads**. Starting read depth recommendations for standard profiling may be adjusted based on user needs.

Archer™ Analysis settings

Sequencing data should be processed using Archer Analysis (v7.0, or greater). The FusionPlex Lung requires selection of the *RNA Fusion* pipeline, found under the *RNA* Input Type (see the Archer Analysis User Guide for more details on setting up your analysis).

Processing of FusionPlex Lung libraries requires a one-time upload of the Panel GTF. Files can be obtained by contacting archer-tech@idtdna.com

Assay targets

Gene	Accession	Exon	Variant Type	Description*
ALK	NM_004304	22, 23, 25	Mutation	T1151-C1156, F1174, L1196-S1206, G1269
ALK	NM_004304	2, 4, 6, 10, 16, 17, 18, 19, 20, 21, 22, 23, 26	Fusion	5'
BRAF	NM_004333	15	Mutation	V600
BRAF	NM_004333	2, 7, 8, 9, 10, 11, 12, 15, 16	Fusion	5'
BRAF	NM_004333	1, 3, 7, 8, 10, 13	Fusion	3'

2425 55th Street, Boulder, CO 80301 | www.archerdx.com | archer-tech@idtdna.com

RA-DOC-028 / REV01

For research use only. Not for use in diagnostic procedures.



Product Insert

FusionPlexTM Lung Panel

0	A •	F	Maniant Ton	D
Gene	Accession	Exon	Variant Type	Description*
EGFR	NM_005228	18, 19, 20, 21	Mutation	E709-G719, E746-L760, V774-G796, L858-L861
EGFR	NM_005228	7, 8, 9, 16, 19, 20,	Fusion	5'
EGFR	NM_005228	8	Exon 2-7 Skipping (EGFRvIII)	5'
EGFR	NM_005228	1, 24, 25	Fusion	3'
EGFR	NM_005228	1	Exon 2-7 Skipping (EGFRvIII)	3'
FGFR1	NM_015850	2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 17	Fusion	5'
FGFR1	NM_015850	12, 17	Fusion	3'
FGFR2	NM_000141	2, 5, 7, 8, 9, 10	Fusion	5'
FGFR2	NM_000141	16, 17	Fusion	3'
FGFR3	NM_000142	3, 5, 8, 9, 10	Fusion	5'
FGFR3	NM_000142	16, 17, 18	Fusion	3'
KRAS	NM_004985	2, 3	Mutation	G12-G13, Q61
MET	NM_000245	2, 4, 5, 6, 13, 14, 15, 16, 17, 21	Fusion	5'
MET	NM_000245	15	Exon 14 Skipping	5'
MET	NM_000245	2, 13	Fusion	3'
MET	NM_000245	13	Exon 14 Skipping	3'
NRG1	NM_013957	1, 8	Fusion	5'
NRG1	NM_004495	1, 2, 3, 4, 6	Fusion	5'
NRG1	NM_013962	1	Fusion	3'

2425 55th Street, Boulder, CO 80301 | www.archerdx.com | archer-tech@idtdna.com

RA-DOC-028 / REV01

For research use only. Not for use in diagnostic procedures.



Product Insert

FusionPlexTM Lung Panel

Gene	Accession	Exon	Variant Type	Description*
NTRK1	NM_002529	2, 4, 6, 8, 10, 11, 12, 13	Fusion	5'
NTRK2	NM_006180	5, 7, 9, 11, 12, 13, 14, 15, 16, 17	Fusion	5'
NTRK3	NM_002530	4, 7, 10, 12, 13, 14, 15, 16	Fusion	5'
NTRK3	NM_001007156	15	Fusion	5'
NTRK3	NM_002530	13, 14, 15	Fusion	3'
RET	NM_020630	15,16	Mutation	A883, M918
RET	NM_020630	2, 4, 6, 8, 9, 10, 11, 12, 13, 14	Fusion	5'
ROS1	NM_002944	38	Mutation	G2032
ROS1	NM_002944	2, 4, 7, 31, 32, 33, 34, 35, 36, 37	Fusion	5'

^{*}The mutations listed under the Description column are targeted by the assay design. Version 6.2 and earlier of Archer Analysis may not support RNA SNV/indel variant calling at exon junctions depending on the sequence context (SNVs ≤5bp, indels ≤30bp). *De Novo* RNA SNV/indel and Internal Tandem Duplication mutation detection are not supported on the Ion Torrent Sequencing System.

2425 55th Street, Boulder, CO 80301 | www.archerdx.com | archer-tech@idtdna.com



Product Insert FusionPlexTM Lung Panel

Limitations of use

For research use only. Not for use in diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

Safety data sheets pertaining to this product are available upon request.

This product or service is licensed under one or more of the following U.S. Patents: 8,835,358; 9,290,808; 9,290,809; 9,315,857; 9,708,659; and 9,816137 owned by BD and is licensed solely for the use described in the associated product literature. No other rights, implied or otherwise, are granted to purchaser hereunder. Purchaser agrees, by way of example and not limitation, not to use this product to trace back the origin of a nucleic acid to an individual cell as a discrete entity (e.g., single cell analysis).

© 2023 Integrated DNA Technologies, Inc. All rights reserved. FusionPlex, VariantPlex, LiquidPlex, Immunoverse, Archer Analysis, and Archer Assay Marketplace are trademarks of Integrated DNA Technologies, Inc. All other marks are the property of their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks.

2425 55th Street, Boulder, CO 80301 | www.archerdx.com | archer-tech@idtdna.com