



Using the rhAmpSeq™ CRISPR Analysis Tool

> SEE WHAT MORE WE CAN DO FOR **YOU** AT **WWW.IDTDNA.COM**.

NGS • CRISPR • Functional Genomics • Synthetic Biology • PCR & qPCR • Custom Oligos • Contract Manufacturing

REVISION HISTORY

Version	Release date	Description of changes
2	June 2023	Updated contact information, account information image, and target report text
1	April 2021	Initial release

Table of contents


Revision history	2
Introduction	4
Access the rhAmpSeq CRISPR Analysis Tool	4
Create an account	4
Manage your account	6
Start analyzing	7
Upload data	7
Upload large files via Connector setup	9
Data with multiple amplicons	10
Quality control checking BED inputs	11
Handling QC statuses	12
Analyzing singleplex experiments	12
Interpreting singleplex experiments	14
Viewing allele details	15
Viewing the target report	17
Screening guides	18
Analyzing multiplex experiments	18
Interpreting multiplex experiments	20
Viewing allele details	21
Viewing the aggregated results	22

INTRODUCTION

This analysis guide describes how to use the rhAmpSeq CRISPR Analysis Tool, which was designed for flexible, cloud-based interrogation of CRISPR-mediated genome editing results. This information can also be accessed in video format by viewing the [rhAmpSeq CRISPR Analysis Tool tutorial](#).


ACCESS THE rhAmpSeq CRISPR ANALYSIS TOOL

From the IDT website, go to Tools, CRISPR Genome Editing, then select the [rhAmpSeq CRISPR Analysis Tool](#).



CREATE AN ACCOUNT

An IDT account provides access to the [Analysis Lab](#) page and the ability to use the rhAmpSeq CRISPR Analysis Tool.



 **Note:** If you do not already have an IDT account, click [Register](#) to create your own account.

Account Information

Login Information

Login Name:

Password:

Confirm Password:

This is the information you will use to log in when you order, view past orders, or change your account information.

Personal Information

First Name:

Last Name:

Email:

Confirm Email:

PI First Name:

PI Last Name:

Institution Type:

Institution Name:

I agree to receive email communications from IDT.:

I am interested in receiving a free consultation call from IDT.: ☐

Important:

Once established, your first and last name cannot be changed without contacting us; for organization or institution changes, we ask that you please create a new account under the new parent organization. Please contact IDT Customer Care if you have any questions or concerns.

IDT uses your contact information in order to provide you with information about our services. IDT does not share your personal information with any non-affiliated third parties, other than our distributors, nor do we process your personal information for any other purpose. To learn more about how IDT treats your information, please review our [privacy statement](#).

SAVE

Once established, you will use this account login to sign in and use the rhAmpSeq CRISPR Analysis Tool.



Search



GET HELP SIGN IN

0 ITEMS \$0.00 USD

PRODUCTS & SERVICES APPLICATIONS & TECHNOLOGIES SUPPORT & EDUCATION TOOLS COMPANY

Please Sign In

UserName

Password

☐ Keep me signed in. [Details](#)

[Forgot Password](#)

SIGN IN

REGISTER

Please sign in to use IDT's custom online ordering tools. If you don't yet have an IDT account, join the IDT community! Create your free account today and enjoy unlimited access to our innovative web tools, streamlined ordering, and expert educational content.

After you log in, the Analysis Lab landing page opens. If you are not already part of an analysis lab, click the **Create a Lab** button and give a unique name to your lab.

Analysis Lab Overview

You don't belong to any active lab. Please join or create one.

CREATE A LAB

REDEEM CREDIT

SEND INVITE

ANALYSIS TOOL

Create New Lab

Lab Name

CRISPAItRationsBFXSupport

CREATE

After creating a lab, you can invite collaborators, colleagues, and students to join your analysis lab by clicking **Send Invite**. Enter a list of those you want to collaborate with, and a custom message to accompany your invitation.

Analysis Lab Overview

CRISPAItRationsBFXSupport

Analysis Credit Balance: 1000

Manage lab

CREATE A LAB

REDEEM CREDIT

SEND INVITE

ANALYSIS TOOL

Members Pending invites

Member Name

Gavin Kurgan

Send Lab Invite(s)

Send to

email address

Invite message

Hi, Join my analysis lab!

SEND



Tip: You can control user privileges by setting permission levels for each invited member.

MANAGE YOUR ACCOUNT

To use the rhAmpSeq CRISPR Analysis Tool, you will need Analysis Credits. You received a code for these credits when the order for your purchased product was invoiced. To redeem your Analysis Credits, click the gray **Redeem Credit** button, then enter the code. View your current balance and transaction history in the top left corner.



Tip: Refresh the page to see updated information on your account balance.

Analysis Lab Overview

CRISPAItRationsBFXSupport

Analysis Credit Balance: 1000

VIEW RECENT TRANSACTIONS

Manage lab

CREATE A LAB

REDEEM CREDIT

SEND INVITE

ANALYSIS TOOL

Members Pending invites

Member Name

Gavin Kurgan

Ellen Black

Redeem Credit


Credit code

IDTcustomer00001


REDEEM

START ANALYZING

Open the rhAmpSeq CRISPR Analysis Tool by clicking the blue **Analysis Tool** button. Select **Confirm** when prompted.



A new login page will open where you will enter your IDT account credentials again for security purposes. Now you are ready to begin analyzing data!

 **Tip:** Since you have learned how to create and manage your account, you are ready to explore the other [tutorial videos](#) that further describe how to use the IDT rhAmpSeq CRISPR Analysis Tool.

At any point if you need help or have questions, please [contact us](#). We are here to assist you.

UPLOAD DATA

Before we delve into the process of uploading data, let us review the supported data formats. Specifically, you will have to provide demultiplexed sequencing data, usually in the form of .fastq or .fastq.gz files as shown in the below image. These files contain the following attributes for each read in the sequenced sample:

- read name
- sequence
- associated quality score



There are a few different ways to upload your data to the rhAmpSeq CRISPR Analysis Tool.


You can:

- drag and drop files into the application
- upload files from a local drive
- upload files from the cloud using the Connector

For smaller file sizes, using drag-and-drop is the fastest and easiest method for uploading sequence data.


Simply navigate to the Files tab, then drag the samples you want analyzed into the open gray box in the interface as shown below.

Files are uploaded in parallel, and upon successful upload, R1/R2 files are automatically paired to create a "Sample."




Notes:

- For drag-and-drop uploads you must stay on the upload page until the transfer is completed. Wait for the progress bars to indicate when the uploads are complete before moving to the next step.
- Alternatively, you can manually pair read data into a sample by selecting the uploaded files, then clicking the orange Create sample button.



UPLOAD LARGE FILES VIA CONNECTOR SETUP


For larger files, or high-throughput applications, the drag-and-drop method may not be fast enough. Consider using the Connector setup upload feature. Navigate to the **Settings** page, indicated by a gear icon.




Here, you can set up the uploading of data from your local drive, or cloud storage systems. This guide will give you a high level overview of how to set up the Connector; for a more in-depth review, see the [rhAmpSeq CRISPR Analysis Tool Data Transfer Guide](#).

To start uploads from your local drive to the analysis platform, click the **New** button under Local Connector Settings, then specify the following:

- name of the connector
- computer's operating system
- upload/download location of your data



Save your inputs, then download and install the connector software.



Follow the on-screen installation prompts.

To start uploads from a cloud streaming scenario, click the **New** button under the appropriate cloud service provider that you use. The setup is the same as above (uploads from a local drive), except that some specifics for each cloud provider may differ. Reference the [rhAmpSeq CRISPR Analysis Tool Data Transfer Guide](#) for detailed instructions based on your cloud provider.

Using the rhAmpSeq CRISPR Analysis Tool

DATA WITH MULTIPLE AMPLICONS

If you are analyzing data with multiple amplicons in each sample, you will need to upload two 6-column BED-formatted files specifying amplicon and guide coordinates of the genome of interest. These BED files contain the following for all targets being analyzed:


- chromosome
- start location
- stop location
- name
- strand

Important: Guide coordinates should not include the protospacer adjacent motif (PAM) but should include the correct strand information (+ or -) so that the guide is 5' to 3'.


Note: If the [rhAmpSeq Design Tool](#) was used to generate the panel, then the amplicon BED file can be found in the Design Tool dashboard, as shown below. The file "Assay_panel.bed" will work for this.


View of GeneX_guides.bed					
File Name		GeneX_guides.bed		Size	3 KB
				Creation Date	2021/01/04 16:49
1	chr11	36573321	36573341	GeneX	0 -
2	chr1	3592654	3592674	OT1	0 -
3	chr1	16415506	16415526	OT2	0 +
4	chr1	33698343	33698363	OT3	0 +
5	chr1	34144141	34144161	OT4	0 +
6	chr1	49123705	49123725	OT5	0 -
7	chr1	52927734	52927754	OT6	0 +
8	chr1	115105040	115105060	OT7	0 -
9	chr1	154428110	154428130	OT8	0 +
10	chr1	166894026	166894046	OT9	0 -
11	chr1	183144152	183144172	OT10	0 -
12	chr1	200374008	200374028	OT11	0 -
13	chr1	206221475	206221495	OT12	0 +
14	chr10	6437025	6437045	OT13	0 -
15	chr10	106920200	106920220	OT14	0 -
16	chr10	121109396	121109416	OT15	0 -
17	chr10	125696309	125696329	OT16	0 +
18	chr11	2525645	2525665	OT17	0 +
19	chr11	12125984	12126004	OT18	0 -
20	chr11	20102285	20102305	OT19	0 +
21	chr11	57348913	57348933	OT20	0 -
22	chr11	64787970	64787990	OT21	0 -
23	chr11	103863224	103863244	OT22	0 -
24	chr12	26976538	26976558	OT23	0 +
25	chr12	49647269	49647289	OT24	0 -
26	chr12	56088641	56088661	OT25	0 +
27	chr12	57198911	57198931	OT26	0 +
28	chr12	57289518	57289538	OT27	0 +
29	chr12	74842086	74842106	OT28	0 +
30	chr12	90948235	90948255	OT29	0 -
31	chr12	106046735	106046755	OT30	0 +
32	chr12	109286629	109286649	OT31	0 -
33	chr12	124489776	124489796	OT32	0 +
34	chr13	22458937	22458957	OT33	0 +


To upload these files, use the gear icon to open the **Settings** page, then scroll down to find the section called "BED Files for Multiplex Analysis." Upload your guide RNA and amplicon BED files into the specified sections by using the drag-and-drop upload method.


IDT
INTEGRATED DNA TECHNOLOGIES


Files

 Samples

 Analyses

 Aggregations

 Settings



<

QUALITY CONTROL CHECKING BED INPUTS

In order to use your BED files for multiplex analysis, the uploaded BED files need to be checked for errors. Navigate to the **Settings** page, then look for the “BED Files Pairs” section. This is where you QC-check your inputs to make sure they can successfully be used for analysis. To begin, select the orange **Pair BED files** button, as shown below.

Settings

⌵ BED Files Pairs

Save

Pair Name	Guide RNA	Amplicons	Species	Creation Date	Status
HBB_guides_dual_min_min	HBB_guides_dual_min_min.bed	HBB_test_amplicons.bed	Homo sapiens - human (GRCh38)	2020/12/04 11:13	Fail
HBB_guides_dual_min_pos	HBB_guides_dual_min_pos.bed	HBB_test_ampliconsFixed.bed	Homo sapiens - human (GRCh38)	2020/12/04 11:57	Pass
HBB_guides_all_dropout_fix	HBB_guides_all_dropout.bed	HBB_amplicons_all_dropout.bed	Homo sapiens - human (GRCh38)	2020/12/04 15:06	Warning
HBB_guides_dual_pos_pos	HBB_guides_dual_pos_pos.bed	HBB_test_ampliconsFixed.bed	Homo sapiens - human (GRCh38)	2020/12/04 14:59	Fail
HBB_guides_dropout_fix	HBB_guides_dropout.bed	HBB_amplicons_dropout.bed	Homo sapiens - human (GRCh38)	2020/12/04 15:08	Warning
HBB_guides_singleplex	HBB_guides_singleplex.bed	HBB_amplicons_singleplex.bed	Homo sapiens - human (GRCh38)	2020/12/04 15:05	Warning
HBB_guides_dual_pos_min	HBB_guides_dual_pos_min.bed	HBB_test_ampliconsFixed.bed	Homo sapiens - human (GRCh38)	2020/12/04 14:58	Pass
HBB_guides_minus	HBB_guides_minus.bed	HBB_test_ampliconsFixed.bed	Homo sapiens - human (GRCh38)	2020/12/04 15:01	Warning
HBB_guides_positives	HBB_guides_positive.bed	HBB_test_ampliconsFixed.bed	Homo sapiens - human (GRCh38)	2020/12/04 15:00	Warning
HBB_guides_dropout	HBB_guides_dropout.bed	HBB_test_ampliconsFixed.bed	Homo sapiens - human (GRCh38)	2020/12/04 15:04	Warning
HBB_guides_multi	HBB_guides_multi.bed	HBB_test_ampliconsFixed.bed	Homo sapiens - human (GRCh38)	2020/12/04 15:03	Warning
HBB_guides_four	HBB_guides_four.bed	HBB_test_ampliconsFixed.bed	Homo sapiens - human (GRCh38)	2020/12/04 14:59	Warning
HBB_guides_dual_min_min_pass	HBB_guides_dual_min_min.bed	HBB_test_ampliconsFixed.bed	Homo sapiens - human (GRCh38)	2020/12/04 11:34	Pass
HBB_guides_all_dropout	HBB_guides_all_dropout.bed	HBB_test_ampliconsFixed.bed	Homo sapiens - human (GRCh38)	2020/12/04 15:04	Fail
DemoPairing	HBB_guides_minus.bed	HBB_amplicons_singleplex.bed	Homo sapiens - human (GRCh38)	2021/02/01 14:34	Ongoing
HBB_guides_positive_1.0.0	HBB_guides_positive.bed	HBB_test_ampliconsFixed.bed	Homo sapiens - human (GRCh38)	2021/02/01 13:57	Warning

Pair BED files

You will enter:

- a unique name to reference the files being paired
- the pipeline version intended for analysis
- guide and amplicon BED files
- species these coordinates come from

These entries will initiate the QC check on the files, which will take approximately 5–10 minutes to complete. A status message of PASS, WARNING, or FAIL will appear when the pairing is complete.

```
02/01/2021 08:20:49 PM - INFO - create:set_up_logging - Logging to bedQC.log
02/01/2021 08:20:49 PM - WARNING - columnBiasCheck:columnBiasCheck - There is only one value of {'+'} in column 5; make sure this .
02/01/2021 08:20:50 PM - INFO - create:set_up_logging - Logging to bedQC.log
```

HANDLING QC STATUSES

If a **PASS** status is shown, you are ready to start your analysis.

If a **WARNING** status is shown, double-click the pairing to understand why the paired BED files were flagged. These files will still be available for analysis; however, they should be reviewed to make sure there are no persistent problems with your experimental inputs.

If a **FAIL** status is shown, there is something wrong with the BED inputs that needs correction. Check the log and make the appropriate corrections in the files to reverse the failure.

Important: Once paired, BED files cannot be deleted from the analysis tool. This is to ensure they are available for previous analyses that used them. Try to provide unique and memorable names for these BED file pairs.

ANALYZING SINGLEPLEX EXPERIMENTS

If you want to screen on-target guides and conditions, you will need to know how to analyze your singleplex experiment.

Important: Make sure that data has been completely uploaded and put into samples before proceeding.

First, select your samples. Navigate to the **Samples** page where you will be able to see all the samples that have FASTQ files associated with them. Select the groups of samples with a single amplicon that you want to analyze. Begin the analysis by clicking **Analyze Singleplex**.

To start the run, you will need to associate information and parameters with the run and individual samples. This is where you will also specify the version of the CRISPRatRations software to use for the analysis.



Before starting the run, provide the information highlighted by red boxes, namely the nucleotide sequences for the guide and the amplicon, 5' to 3', excluding the PAM.

Sample	Guide Sequence	Amplicon Sequence	Nuclease	Window Size	Donor
Singleplex_GeneX_Control	<input type="text"/>	<input type="text"/>	Cas9	default	<input type="text"/>
Singleplex_GeneX_EditedHDR	<input type="text"/>	<input type="text"/>	Cas9	default	<input type="text"/>
Singleplex_GeneX_EditedOnly	<input type="text"/>	<input type="text"/>	Cas9	default	<input type="text"/>

The additional parameters that can be configured on a run include:

- the applicable CRISPR nuclease
- window size for interrogating variants
- nucleotide sequences for HDR donors
- well position (for advanced visualizations)
- sample name
- treated vs. untreated sample

Tip: You can quickly populate the same entry down an entire column at once by entering the information in the first row, then clicking the down arrow icon to the right of the field.

Note: For batch runs that have a large number of samples, you will need to upload a samplesheet to provide run information. To get a template of this samplesheet, click the **Export to CSV** button.

Sample	Guide Sequence	Amplicon	Nuclease	Window	S Donor	Well	Group ID	Treated
Singleplex_GeneX_Control	GGCATCGTGTACTACACGG	GCTCAGA: Cas9	default			A1	Control	FALSE
Singleplex_GeneX_EditedHDR	GGCATCGTGTACTACACGG	GCTCAGA: Cas9	default		ATGGTCT	A3	EditedHDR	TRUE
Singleplex_GeneX_EditedOnly	GGCATCGTGTACTACACGG	GCTCAGA: Cas9	default			A2	Edited	TRUE

After specifying your sample information in Microsoft Excel, save and upload this file directly into the interface using the **Import from CSV** button.

Before starting your run, designate an "Aggregation" name, which creates a single location to interrogate all samples in the run. You can add additional metadata and unique descriptions to help label samples and experimental details for tracking results before starting the run.

Caution: When a red box highlights any of the sample information fields, you will not be able to start the run. Check the hover text on any parameter with a red box to find acceptable inputs and make corrections accordingly. When all inputs are satisfied and no red boxes appear, then you are ready to begin the run.

Run Multiplex Analysis

Selected samples (3)
AAVS1_Ctrl, AAVS1_HiFiCas9 and AAVS1_WTCas9

Analysis mode: CRISPAIt-Rations 3.2.4

Analysis Tags: Type to add a tag

Parameters

Sample	BED Files Pair	Nuclease	Window Size	Donor
AAVS1_Ctrl	[Red Box]	Cas9	default	
AAVS1_HiFiCas9		Cas9	default	
AAVS1_WTCas9		Cas9	default	


Delimiter: . [Export to CSV] [Import from CSV]

Tooltip: BED Files Pair is required


To check the status or results of a run, navigate to the **Analyses** page. As soon as cloud resources begin processing, a progress bar should appear. Analyses can take between 20 min–2 hr to complete.

INTERPRETING SINGLEPLEX EXPERIMENTS

View your results by double-clicking the sample on the **Analyses** page. This opens the landing page where you will interpret the results of your experiment. View the graphics for all targets amplified in the **Target Details** section. You can also easily export the graphics in the interface for selected targets using the **Export** button in the top right corner of the page.




Also, on this page you will be able to find the frequency of the editing, the specific repair pathways used, and any frameshifting events that may have occurred. Dive deeper to explore the exact position of insertion and deletion events in the sequenced sample, even the size of these events. We include a frequency histogram showing the normalized positional frequency of insertions, as shown below.



VIEWING ALLELE DETAILS

Go to the **Alignment Details** page to view the distribution of alleles in the sequenced population using the Integrative Genomics Viewer (IGV).

Click on the event of interest to see both the frequency of the event and the read count.



For those with advanced graphical considerations, the raw data we use to render the graphics can be downloaded from the Analysis Details page.

Analysis Details

Target Report

Target Details

Alignment Details

Guide Sequence

GGCATCGTGTAACACGG

Nuclease

Cas9

Max Overlap

152

Min Read Size

60

Window Size

default

Group ID

Edited

Treated

TRUE

Species

notSpecified

Total on-target reads

5,078

% Mapped

99.84

Reads (total)

5170.0

% QC passing

99.45841

% Merged

98.91093

% Primer-dimer

0.0

% On-target

99.84

Results

File Name	Size	Format	Download
OUTPUT.tar.gz	821.27 KB	UNKNOWN	
SAMPLE_allSnpsTables.csv	11.74 KB	CSV	
SAMPLE_AITargetEditing.png	11.69 KB	PNG	
SAMPLE_collapsedIndelInfo.csv	305.19 KB	CSV	
SAMPLE_deleteLocationHist.csv	8.14 KB	CSV	
SAMPLE_indel.csv	2.96 KB	CSV	
SAMPLE_insertLocationHist.csv	5.62 KB	CSV	
SAMPLE_summary.csv	704 B	CSV	
SUMMARIES.tar.gz	17.94 KB	UNKNOWN	

Results files of the targets are consultable in the Target Details menu.

On this page, you can see general sample information, as well as an assortment of other files for either viewing or downloading.

Analysis Details: 20191204-CP3P7-S697, CRISPAIt-Rations 3.2.4, 2021/01/08 09:45

total on-target reads: 12775

View of SAMPLE_collapsedIndelInfo.csv

File Name: SAMPLE_collapsedIndelInfo.csv

Size: 23.99 KB

Creation Date: 2021/01/08 09:55

freqCount	ReferenceAlignment
211	ACAGGACACTTCCTTGCAGGAACAGAGGTGCCATGCAGCCCCGGGTACTCCTTGTGTGTC
109	ACAGGACACTTCCTTGCAGGAACAGAGGTGCCATGCAGCCCCGGGTACTCCTTGTGTGTC
69	ACAGGACACTTCCTTGCAGGAACAGAGGTGCCATGCAGCCCCGGGTACTCCTTGTGTGTC
69	ACAGGACACTTCCTTGCAGGAACAGAGGTGCCATGCAGCCCCGG-- GTACTCCTTGTGTGTCCTCTGGCGCTCTGGCCTCTGCCCCGTAAGCACTTGGTGG
44	ACAGGACACTTCCTTGCAGGAACAGAGGTGCCATGCAGCCCCGGGTACTCCTTGTGTGTC
43	ACAGGACACTTCCTTGCAGGAACAGAGGTGCCATGCAGCCCCGG--


Close

VIEWING THE TARGET REPORT


To interpret and compare results of different samples within an experiment, go to the **Aggregations** tab.

On the **Target Report** landing page, view all the samples that were aggregated for summarization. Any report that you view can easily be exported to Microsoft Excel. From the **Target Report** page, sort or filter the experimental results as needed, and create additional graphics in the software of your choice.

To compare a sample with another, click **Target Details**. On this page you will find the same graphics that were available for individual samples; however, here you can also overlay results from other relevant samples within the whole experiment.



In the graphic below, the run contained two samples with different guides, treated under the same conditions. Notice that the indel profiles of the two guides are different, as would be expected.




SCREENING GUIDES

When screening a large number of guides under varying conditions and treatments, you can use the rhAmpSeq Analysis Tool to generate graphics in 96- or 384-well plate formats to find conditions that are optimal for your experimental goals. The graphics in 96- or 384-well plate formats are available for singleplex analyses only. To create these graphics, simply specify the “Well” information for each singleplex sample prior to sample analysis.

Once your run is complete, go to the **Aggregations** tab where you can find graphics displaying the following:

- frequency of editing
- repair pathway utilization
- frameshift frequency



ANALYZING MULTIPLEX EXPERIMENTS

The most common application for multiplex CRISPR NGS is to quantify on- and off-target editing effects of a guide. In our example, we show how to analyze samples that have multiple amplicons per sample, known as multiplexed samples.

Important: Make sure the following have been completed before proceeding.

- Data have been uploaded and created into samples
- Guide and amplicon BED files have been uploaded and paired

Tip: For a deep dive into pairing your guide and amplicon BED files, watch the section of the rhAmpSeq CRISPR Analysis Tool video [How to quality control check multiplex file inputs](#).

To begin, select the samples you want to analyze, then click **Analyze Multiplex**.

The screenshot shows the 'Samples' tab in the IDT software interface. At the top, there are tabs for 'Files', 'Samples', 'Analyses', and 'Aggregations'. Below the tabs, there are buttons for 'Analyze Singleplex', 'Analyze Multiplex' (which is highlighted with a mouse cursor), 'Edit', 'Unlink Files', and 'Delete'. A table below shows a list of samples with columns for 'Sample Name', 'Creation Date', 'Status', 'Size', and 'Experiment'. The samples listed are AAVS1_Ctrl, AAVS1_WTCas9, and AAVS1_HiFiCas9, all with a status of 'New'.

Sample Name	Creation Date	Status	Size	Experiment
AAVS1_Ctrl	2021/02/02 15:36	New	2.74 MB	Multiplex
AAVS1_WTCas9	2021/02/02 15:35	New	3.13 MB	Multiplex
AAVS1_HiFiCas9	2021/02/02 15:35	New	2.18 MB	Multiplex

Once your data has successfully been checked for quality control, select your named guide/amplicon pairs and other relevant run information, including:

- the appropriate CRISPR nuclease
- window size for interrogating variants
- nucleotide sequences for HDR donors
- sample name information
- treated vs. untreated

Navigate to the **Analyses** page to check the status or results of your analysis. As soon as cloud resources begin processing, you will see the progress bar. Analyses can take anywhere from 20 min–2 hr to complete. Your turnaround time will depend on your read depth plus the number of amplicon targets being quantified per sample.

The screenshot shows the 'Analyses' tab in the IDT software interface. On the left, there is a sidebar with the IDT logo and a section titled 'Available analyses' showing 'CRISPRaltRatios' with a count of 50. The main area displays a table of analysis results with columns for 'Sample', 'Progress', 'Status', 'Reads (total)', '% QC passing', and '% Merged'. The table lists 18 samples, all with a status of 'success'.

Sample	Progress	Status	Reads (total)	% QC passing	% Merged
20191204-CP3P7-S978	██████████	success	4,561	99.76108	99.629875
20191204-CP3P7-S923	██████████	success	33,598	99.65163	99.54672
20191204-CP3P7-S914	██████████	success	10,308	99.50957	99.623116
20191204-CP3P7-S71	██████████	success	6,047	99.523575	99.83493
20191204-CP3P7-S80	██████████	success	21,694	99.619705	99.79303
20191204-CP3P7-S1259	██████████	success	5,841	99.81324	99.64279
20191204-CP3P7-S135	██████████	success	3,389	99.735374	99.97052
20191204-CP3P7-S416	██████████	success	3,398	99.82619	99.44863
20191204-CP3P7-S352	██████████	success	8,849	99.54041	99.71847
20191204-CP3P7-S361	██████████	success	32,353	99.60705	99.71336
20191204-CP3P7-S642	██████████	success	38,877	99.583984	99.64375
20191204-CP3P7-S1204	██████████	success	39,934	99.646835	99.67803
20191204-CP3P7-S697	██████████	success	2,275	99.82832	98.79622
20191204-CP3P7-S1195	██████████	success	4,938	99.53704	99.858444
20191204-CP3P7-S633	██████████	success	7,462	99.469284	99.573166


INTERPRETING MULTIPLEX EXPERIMENTS

When you have completed your analysis, the **Target Details** landing page displays. This page contains graphics for all amplified targets. Each target can be toggled on or off by checking or unchecking the checkboxes to the left of the targets in the table. Also, the active display on the page can easily be exported by using the **Export** button in the top right corner.

From the **Target Details** page, you can find a range of data, including:


- frequency of editing
- repair pathways
- frameshifting events

Additionally, you can explore the position of any insertion and deletion event that occurred in the sequenced sample and the size of these events. We have included a frequency histogram that shows the normalized positional frequency of insertions below.




VIEWING ALLELE DETAILS

Go to the **Alignment Details** page to view the distribution of alleles in the sequenced population using IGV. Click on the event of interest to see both the frequency of the event and the read count.




For those with advanced graphical considerations, the raw data we use to render these graphics can be downloaded from the **Analysis Details** page. On this page, you can see general sample information, as well as an assortment of other files for either viewing or downloading.




AutoSave: Off Target Report: Multiplex 2021-02-02 15:48:21 - Excel

A quick sort by **Target**, **Group ID**, and **Treated** can quickly pair treatments with their corresponding controls to visualize editing differences using the conditionally formatted output.



To investigate differences between samples further, go to the **Target Details** page. Here, you can overlay the results of targets from different experimental samples, then export these graphics using the **Export** button in the top right corner. The graphics provided here are the same as when analyzing a single sample, but single samples cannot be overlaid. An example of this is shown below.



A complete video tutorial version of this guide can be viewed on the [rhAmpSeq CRISPR Analysis Tool](#) page.

Using the rhAmpSeq CRISPR Analysis Tool

For more information, go to: www.idtdna.com/ContactUs

For more than 30 years, IDT's innovative tools and solutions for genomics applications have been driving advances that inspire scientists to dream big and achieve their next breakthroughs. IDT develops, manufactures, and markets nucleic acid products that support the life sciences industry in the areas of academic and commercial research, agriculture, medical diagnostics, and pharmaceutical development. We have a global reach with personalized customer service.

> SEE WHAT MORE WE CAN DO FOR YOU AT WWW.IDTDNA.COM.

For Research Use Only. Not for use in diagnostic procedures. Unless otherwise agreed to in writing, IDT does not intend these products to be used in clinical applications and does not warrant their fitness or suitability for any clinical diagnostic use. Purchaser is solely responsible for all decisions regarding the use of these products and any associated regulatory or legal obligations.

© 2023 Integrated DNA Technologies, Inc. All rights reserved. Trademarks contained herein are the property of Integrated DNA Technologies, Inc. or their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks.
Doc ID: RUO23-1915_001 06/23