

Delivery of ribonucleoprotein complexes into *P. falciparum*-infected erythrocytes by electroporation

Contributed by Dr Simon Cobbold, Department of Biochemistry and Molecular Biology, Bio21 Institute, University of Melbourne, Australia

The method presented here is provided by customers who have used the Alt-R CRISPR-Cas9 System. This method can serve as a starting point for using the Alt-R CRISPR-Cas9 System in similar model organisms but may not be fully optimized for your gene or application. IDT does not guarantee these methods, and application specialists at IDT can only provide general guidance with limited troubleshooting and support.

Consumables and equipment

Consumables from IDT

Item	Catalog #
Alt-R CRISPR-Cas9 crRNA*	IDT predesigned and custom crRNA*: www.idtdna.com/CRISPR-Cas9
Alt-R CRISPR-Cas9 tracrRNA	1072532, 1072533, 1072534
Alt-R S.p. Cas9 Nuclease V3†	1081058, 1081059
Alternatives:	
Alt-R S.p. HiFi Cas9 Nuclease V3	1081060, 1081061
Alt-R S.p. Cas9 D10A Nickase V3	1081062, 1081063
Alt-R S.p. Cas9 H840A Nickase V3	1081064, 1081065
Nuclease-Free IDTE Buffer, pH 7.5 (1X TE solution)	11-01-02-02

* We guarantee the performance of our predesigned gRNAs targeting human, mouse, rat, zebrafish, or nematode genes. For other species, use our proprietary algorithms to design custom gRNAs. If you have spacer designs of your own or from publications, use our design checker tool to assess their on- and off-targeting potential before ordering gRNAs that are synthesized using our Alt-R gRNA modifications. For details about the predesigned gRNA guarantee, see www.idtdna.com/CRISPR-Cas9.

† Alt-R S.p. Cas9 Nuclease V3 (wild-type) is suitable for most genome editing studies. However, some experiments may benefit from use of Alt-R S.p. HiFi Cas9 Nuclease V3, which has been engineered to reduce off-target effects, while retaining the on-target potency of wild type Cas9. Alt-R Cas9 nickases create single-stranded breaks. When a nickase variant is used with 2 gRNAs, off-target effects are reduced, and homology-directed repair is promoted.

Consumables from other suppliers

Item	Supplier	Catalog #
Absolute ethanol	General laboratory supplier	Varies
Cytomix transfection buffer	Custom: 120 mM KCl, 0.15 mM CaCl ₂ , 2 mM EGTA, 5 mM MgCl ₂ , 10 mM K ₂ HPO ₄ /KH ₂ PO ₄ , pH 7.6, 25 mM HEPES, pH 7.6	N/A
Delivery/Rescue plasmid	Custom (context-dependent)	N/A
Fresh O+ Blood	Red Cross / Donation	N/A
RPMI 1640 GlutaMAX™ media	Thermo Fisher Scientific	72400-047
10% AlbuMAX™ II Lipid-Rich BSA (w/v)	Thermo Fisher Scientific	11021045
Hypoxanthine	Sigma/Merck	H9377
Gentamycin sulfate	Sigma/Merck	G1914
Sodium acetate, 3M	Sigma/Merck	S2889

Equipment

Item	Supplier	Catalog #
Gene Pulser Xcell™ Electroporation System	Bio-Rad	1652660
Electroporation Cuvette Plus (400 µL – 2 mm)	BTX	45-0125

Protocol

Culture *P. falciparum*-infected erythrocytes

1. Maintain *P. falciparum*-infected erythrocyte cultures under standard culturing conditions and routinely synchronize with 5% sorbitol.
2. Perform RNP transfections on ring-stage cultures at 6–10% parasitemia (4% haematocrit).

Prepare RNA

1. Resuspend the crRNA and tracrRNA in Nuclease-Free IDTE Buffer:
 - Alt-R CRISPR-Cas9 crRNA, 2 nmol in 20 μL (100 μM)
 - Alt-R CRISPR-Cas9 tracrRNA, 5 nmol in 50 μL (100 μM)
2. Mix the crRNA and tracrRNA oligos in equimolar concentrations (1 μL each, per transfection) in a sterile microcentrifuge tube to a final duplex concentration of 50 μM .
3. Incubate at 95°C for 5 minutes and then let cool to room temperature on the bench.

Form RNP complex

1. Mix 1.4 μL of the RNA solution with 1 μL of Cas9 enzyme per transfection (provided as a 61 μM suspension, 10 $\mu\text{g}/\mu\text{L}$).
2. Incubate at room temperature for 20 minutes.

Prepare rescue plasmid/RNP complex

1. Linearize 50 μg of user-provided plasmid (per transfection) using two compatible restriction enzymes and incubate overnight.



Note: The rescue plasmid is dependent on your desired modification; however, linearizing the plasmid prevents the need for drug cycling or other methods for plasmid elimination following transfection/parasite recovery.

2. Precipitate the linearized rescue template by adding 3M sodium acetate to a final concentration of 10% (v/v), then add ethanol to a final concentration of 70% ethanol (v/v).
3. Incubate at -20°C for at least 4 hours.
4. Centrifuge at 14,000 $\times g$ for 15 minutes (at 4°C).
5. Wash the DNA pellet with ice-cold 70% ethanol and centrifuge again.
6. Using a fresh pipette tip, remove the supernatant.

7. Dry the linearized rescue template for 3 minutes in a sterile environment, then resuspend in 400 μL of cytomix.
8. Mix 2 μL (up to 8 μL) of RNP complex with the rescue template suspension.



Tip: Typically, 2 μL is sufficient to recover integrated parasites; however, you may need up to 8 μL of RNP complex for a standard transfection.

Transfect cells

1. Prepare complete RPMI 1640 media by adding 10% AlbuMAX II (w/v) to a final concentration of 0.5%, hypoxanthine to 200 μM , and gentamycin sulfate to 20 mg/L.
 2. Transfer 1.5 mL of ring-stage *P. falciparum*-infected erythrocyte culture (4% haematocrit, 6–10% parasitemia) to a 10 mL tube.
 3. Centrifuge the cell suspension for 5 minutes at 500 x *g*.
 4. Aspirate the supernatant and resuspend the cell pellet with the RNP-Rescue plasmid solution (400 μL).
 5. Transfer the cell suspension to an electroporation cuvette.
 6. Start electroporation with the following parameters:
 - Voltage: 310 V
 - Capacitance: 950 μF
 - Resistance: ∞
-  **Note:** Run electroporation for approximately 8–12 msec.
7. Immediately transfer the cell suspension to a fresh 10 mL tube with complete RPMI 1640.
 8. Centrifuge the cell suspensions for 5 minutes at 500 x *g*.
 9. Using a fresh pipette tip, aspirate the supernatant.
 10. Add 10 mL of RPMI and 200 μL of fresh O+ blood.
 11. Store the cultures overnight under standard culturing conditions.
 12. Apply the appropriate drug selection the following day.
 13. Monitor the cultures by Giemsa smear, then replace the media every other day until *P. falciparum* is detectable (typically 2–3 weeks).

Because *P. falciparum* does not contain NHEJ, no indels form at the PAM site, thus explains why you will always have to provide a rescue template for DNA recombination. Now you are ready to perform a PCR to confirm integration.

Delivery of ribonucleoprotein complexes
into *P. falciparum*-infected erythrocytes
by electroporation

Integrated DNA Technologies, Inc. (IDT) is your Advocate for the Genomics Age. For more than 30 years, IDT's innovative tools and solutions for genomics applications have been driving advances that inspire scientists to dream big and achieve their next breakthroughs. IDT develops, manufactures, and markets nucleic acid products that support the life sciences industry in the areas of academic and commercial research, agriculture, medical diagnostics, and pharmaceutical development. We have a global reach with personalized customer service. See what more we can do for **you** at www.idtdna.com.

Technical support:
applicationsupport@idtdna.com

For Research Use Only.

© 2019 Integrated DNA Technologies, Inc. All rights reserved. Alt-R is a trademark of Integrated DNA Technologies, Inc. and registered in the USA. All other trademarks are the property of their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks. CRS-10155-PR 08/2019