

Product: xGen™ CDI Primers, 96 rxn

Product number: 10009815

Batch number: 0000777004

Expiration date: 2024-NOV-22

Analysis date: 2021-MAR-31

Component	Component Part Number	Component Batch Number
xGen CDI Primers Index D701, 8 rxn	10012218	0000773275
xGen CDI Primers Index D702, 8 rxn	10012219	0000773276
xGen CDI Primers Index D703, 8 rxn	10012220	0000773278
xGen CDI Primers Index D704, 8 rxn	10012221	0000773279
xGen CDI Primers Index D705, 8 rxn	10012222	0000773280
xGen CDI Primers Index D706, 8 rxn	10012223	0000773281
xGen CDI Primers Index D707, 8 rxn	10012224	0000773282
xGen CDI Primers Index D708, 8 rxn	10012225	0000773283
xGen CDI Primers Index D709, 8 rxn	10012226	0000773284
xGen CDI Primers Index D710, 8 rxn	10012227	0000773285
xGen CDI Primers Index D711, 8 rxn	10012228	0000773286
xGen CDI Primers Index D712, 8 rxn	10012229	0000773287
xGen CDI Primers Index D501, 12 rxn	10012210	0000773267
xGen CDI Primers Index D502, 12 rxn	10012211	0000773268
xGen CDI Primers Index D503, 12 rxn	10012212	0000773269
xGen CDI Primers Index D504, 12 rxn	10012213	0000773270
xGen CDI Primers Index D505, 12 rxn	10012214	0000773271
xGen CDI Primers Index D506, 12 rxn	10012215	0000773272
xGen CDI Primers Index D507, 12 rxn	10012216	0000773273
xGen CDI Primers Index D508, 12 rxn	10012217	0000773274

Parameter	Specification	Result
Functional Testing	Each set of reagents is functionally tested with xGen DNA Library Prep EZ kit and compared to established performance metrics. Libraries are prepared using high quality human gDNA fragmented to 200 bp. Libraries are sequenced on an Illumina® sequencer. Metrics evaluated include library yields, pass filter clusters, average percent Q30 scores, reads in the undetermined bin.	Pass

Storage: Store xGen CDI Primers in a sealed container at –20°C.

Verified by: Krystal Mullnix

Quality release date: 2022-DEC-01

IDT verifies that the information contained herein is true and correct to the best of our knowledge. This document was produced electronically and is valid without signature.

IDT is a trademark of Integrated DNA Technologies, Inc.
Illumina is a registered trademark of Illumina, Corp.