

Archer™ MBC Adapters A9 - A16 for Illumina® Platform

PI-SA0041

Product Description

The Archer Molecular Barcode (MBC) Adapters are used in conjunction with the Archer FusionPlex™ or VariantPlex™ kits in order to create target-enriched libraries for next-generation sequencing on the MiSeq® and NextSeq® platforms.

Archer Illumina assays utilize dual-index barcoding to distinguish between samples. The first index is added just before Adapter Ligation and is embedded in the Archer MBC Adapters. The second index is added in Second PCR and is embedded in the Illumina Index 1 Primers within the Second PCR reaction.

In addition to this index, Archer MBC Adapters contain molecular barcodes – unique nucleotide sequences that can be used to differentiate individual nucleic acid molecules in a sample. Using Archer Analysis, individual MBCs in a sequencing reaction can be identified and the reads de-duplicated. MBCs can also be used to correct for late-cycle PCR artifacts or sequencing errors, enabling high-confidence calling of gene fusions, SNPs and other mutations.

For Research Use Only. Not for use in diagnostic procedures.

Kit Contents

Part #	Name	Quantity	Storage Conditions
SA0041	Barcode Adapters A9 - A16 for Illumina Platform	Barcoded adapters with a total of 8 different index sequences	2°C to 8°C

Additional Information

	MBC Adapter	Index
Illumina Barcode Adapters A9 - A16 (SA0041)	Illumina_RBC1_BC9	GAGGCTTA
	Illumina_RBC1_BC10	GACTTCGT
	Illumina_RBC1_BC11	GTCGTTGT
	Illumina_RBC1_BC12	TCTGCTTC
	Illumina_RBC1_BC13	GAGGACTT
	Illumina_RBC1_BC14	TCCTCGAT
	Illumina_RBC1_BC15	AAGCAGGA
	Illumina_RBC1_BC16	TTCTAGGC



Limitations of Use

For Research Use Only. Not for use in diagnostic procedures.

This product was developed, manufactured, and sold for in vitro use only. The product is not suitable for administration to humans or animals. SDS sheets relevant to this product are available upon request.

© 2016 ArcherDX, Inc. All rights reserved. Archer™, FusionPlex™, and VariantPlex™ are trademarks of ArcherDX, Inc. Illumina®, NextSeq® and MiSeq® are registered trademarks of Illumina, Inc.

For more information please visit <http://www.archerdx.com>



ArcherDX, Inc.
2477 55th Street, Suite 202
Boulder, CO 80301
303-357-9001
www.archerdx.com